

**FORENSIC TOXICOLOGY LABORATORY  
OFFICE OF CHIEF MEDICAL EXAMINER  
NEW YORK CITY**

**HYDROLYSIS OF SPECIMENS FOR SOLID PHASE EXTRACTION**

**PRINCIPLE**

Many drugs are metabolized by the process of glucuronidation. Glucuronides do not extract well in most common analytical procedures. This procedure is designed to hydrolyze the bond between a drug and glucuronic acid resulting in an increase in the amount of unconjugated (free) drug available for extraction from biological specimens. A single point calibrator at a concentration of 1000 ng/mL, a control at 200 ng/mL, a matrix blank and a hydrolysis control consisting of morphine glucuronide spiked into negative matrix to a concentration of 80 ng/mL, is run with each batch. Refer to the Opiate-BE procedure for preparation of additional solutions and standards as required.

**SAFETY**

The handling of all reagents, samples and equipment is performed within the guidelines which are detailed in the safety manual.

**REAGENTS AND MATERIALS**

1. **Deionized water**
2. **HCl, concentrated, 12.1 N** Fisher Scientific or equivalent
3. **100 mM phosphate buffer (pH 6.0)**

Dissolve 3.4 g  $\text{Na}_2\text{HPO}_4$  and 24.2 g  $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$  in 1600 mL DI  $\text{H}_2\text{O}$ .

Dilute to 2000 mL using DI  $\text{H}_2\text{O}$ . Mix. Adjust pH to 6.0 + 0.1 with 100 mM monobasic sodium phosphate (lowers pH) or 100 mM dibasic sodium phosphate (raises pH).

Store at 5 °C in glass.

Stability: 1 month. Inspect each day before use for contamination.

4. **6N HCl**

Add concentrated HCl to water (1:1). Mix thoroughly.

**CAUTION: Prepare this solution in a fume hood. Use appropriate safety equipment.**

5. **Hydrolysis control solution**

Morphine glucuronide (10 mg/L)

Pipet 1 mL of 1.0 mg/mL morphine glucuronide into a 100 mL volumetric flask. Q.S. to 100 mL with methanol. Transfer to properly labeled container.

## PROCEDURE

1. In properly labeled 16 x 150 screw-cap centrifuge tubes, place 1 mL samples, calibrator, and controls as described above, fortified with internal standard. Add a hydrolysis control by pipetting 1.0 ml of blank matrix to a properly labeled tube and fortifying with 40  $\mu$ L of 10 mg/L morphine glucuronide control solution and the appropriate amount of internal standard solution.
2. Add 0.5 mL 6N HCl to all tubes. Vortex to mix Place the tube in a boiling water bath for 45 minutes. Remove the tube and allow it to cool to room temperature.
3. Add 2.0 mL deionized H<sub>2</sub>O and 2.0 mL 100 mM phosphate buffer pH 6.0, cap, and mix by Vortex for 30 seconds. Centrifuge for 15 minutes at  $\approx$  3000 RPM
4. Filter sample through Whatman filter paper into a labeled 25 mL beaker.
5. Decant the supernatant into appropriately labeled Polychrom Clin II column and apply nitrogen at a pressure of 2-4 psi.
6. Proceed with analysis as in Opiate-BE procedure.