

# **Forensic Molecular Biology**

## **Solutions Manual**

**Version 4.0**

Initials: *RCJ*

Date: *2/9/98*

# Alphabetical Table of Contents

S113 Acrylamide Solution, 40%	6
S115 Ammonium Persulfate (0.5g Aliquot)	7
S018 Analytical Gel Loading Buffer	8
S131 Blue Formamide + Loading Buffer (5:1)	8A
S117 BSA Solution, 5 mg/mL	9
S060 Calibration Control	10
S010 Cell Lysis Buffer (CLB)	12
S064 Cell Pellet Control	13
S022 Chelex, 5%	15
S082 Chelex, 20%	17
S027 Chloroform-Isoamyl Alcohol	18
S104 Chromogen Solution	19
S013 Denaturation Solution	20
S094 Digest Buffer	21
S061 Digestion Control	22
S003 DQ $\alpha$ Citrate Buffer	24
S004 DQ $\alpha$ Hybridization Solution	25
S005 DQ $\alpha$ Wash Solution	26
S093 DTT, 1M	27
S030 DTT, 0.39M	28
S009 EDTA, 0.5M	29
S110 EDTA, 200 mM	30
S056 Ethanol, 70%	31
S031 Gel Neutralization Buffer	32
S116 GS500 + Loading Buffer (1:3)	33
S016 Hae III Buffer, 10X	34
S102 Heparin Solution	35
S026 Hybridization Solution, RFLP	36
S079 Hydrogen Peroxide, 3%	37
S032 Lambda Marker	38
S015 Lithium Chloride, 7.5 M	40
S055 Magnesium Chloride, 0.1M	41
S112 dNTP's (10mM stock)	42
S006 Phenol	43
S053 Phenol-Chloroform-Isoamyl Alcohol	45
S042 Phi-X Marker	46
S034 Phosphate Buffered Saline (PBS)	48
S109 PM PCR Reaction Mixture	49
S126 Positive Control-External	50
S127 Positive Control-Internal	51
S097 Pre-Wetting Solution	53
S118 Primer, VWA 1 (50 $\mu$ M)	54
S119 Primer, VWA 2 (50 $\mu$ M)	56
S120 Primer THO1 1 (50 $\mu$ M)	58
S121 Primer THO1 2 (50 $\mu$ M)	60
S122 Primer F13A1 1 (50 $\mu$ M)	62

Initials: *RCJ*Date: *2/9/98*

S123 Primer, F13A1 2 (50 $\mu$ M) .....	64
S124 Primer, FES/FPS 1 (50 $\mu$ M) .....	66
S125 Primer FES/FPS 2 (50 $\mu$ M) .....	68
S011 Protein Lysis Buffer (PLB) .....	70
S014 Proteinase-K Enzyme, 10mg/mL .....	71
S037 Proteinase-K Enzyme, 20mg/mL .....	72
S111 QUAD STR/PCR Reaction Mixture .....	73
S100 QuantiBlot DNA Standards .....	75
S099 QuantiBlot Wash Solution .....	77
S054 Sarkosyl, 10% .....	78
S040 Sarkosyl, 20% .....	79
S132 Sequencing Loading Buffer .....	79A
S101 SDS, 0.1% .....	80
S045 SDS, 10% .....	81
S001 SDS, 20% .....	82
S046 SLS, 20% .....	83
S043 Sodium Acetate, 2 M .....	84
S044 Sodium Acetate, 0.2M .....	85
S035 SP, 25X .....	86
S108 SP, 2X .....	87
S098 Spotting Solution .....	88
S002 SSPE, 20X .....	89
S047 Stain Extraction Buffer .....	90
S059 Sterile Water .....	91
S023 Stripping Solution .....	92
S017 TRIS-HCl, 1M - PH 8.0 .....	93
S128 TE <sup>-4</sup> , 1X .....	94
S050 Test Gel Loading Buffer .....	95
S062 Test Gel Standard .....	96
S051 TNE, 10X .....	98
S107 TNE, 1X .....	99
S052 TRIS, 0.1M - pH 7.8 .....	100
S114 Urea (25 g Aliquot) .....	101
S114a (18 g Aliquot-377 Sequencer) .....	101A
S114b Urea Urea (10.8 g Aliquot-377 Sequencer) .....	101B
S024 Wash Solution #1 .....	102
S025 Wash Solution #2 .....	103
S020 Yield Calibrators .....	104
S021 Yield Gel Loading Buffer .....	106

Initials: *RCJ*

Date: *2/9/98*

## Numerical Table of Contents

S001 SDS, 20%	82
S002 SSPE, 20X	89
S003 DQ $\alpha$ Citrate Buffer	24
S004 DQ $\alpha$ Hybridization Solution	25
S005 DQ $\alpha$ Wash Solution	26
S006 Phenol	43
S009 EDTA, 0.5M	29
S010 Cell Lysis Buffer (CLB)	12
S011 Protein Lysis Buffer (PLB)	70
S013 Denaturation Solution	20
S014 Proteinase-K Enzyme, 10mg/mL	71
S015 Lithium Chloride, 7.5 M	40
S016 Hae III Buffer, 10X	34
S017 TRIS-HCl, 1M - PH 8.0	95
S018 Analytical Gel Loading Buffer	8
S020 Yield Calibrators	104
S021 Yield Gel Loading Buffer	106
S022 Chelex, 5%	15
S023 Stripping Solution	92
S024 Wash Solution #1	102
S025 Wash Solution #2	103
S026 Hybridization Solution, RFLP	36
S027 Chloroform-Isoamyl Alcohol	18
S030 DTT, 0.39M	28
S031 Gel Neutralization Buffer	32
S032 Lambda Marker	38
S034 Phosphate Buffered Saline (PBS)	48
S035 SP, 25X	86
S037 Proteinase-K Enzyme, 20mg/mL	72
S040 Sarkosyl, 20%	79
S042 Phi-X Marker	46
S043 Sodium Acetate, 2 M	84
S044 Sodium Acetate, 0.2M	85
S045 SDS, 10%	81
S046 SLS, 20%	83
S047 Stain Extraction Buffer	90
S050 Test Gel Loading Buffer	95
S051 TNE, 10X	98
S052 TRIS, 0.1M - pH 7.8	100
S053 Phenol-Chloroform-Isoamyl Alcohol	45
S054 Sarkosyl, 10%	78
S055 Magnesium Chloride, 0.1M	41
S056 Ethanol, 70%	31
S059 Sterile Water	91
S060 Calibration Control	10
S061 Digestion Control	22

Initials: *RG*Date: *2/9/98*

S062 Test Gel Standard	96
S064 Cell Pellet Control	13
S079 Hydrogen Peroxide, 3%	37
S082 Chelex, 20%	17
S093 DTT, 1M	27
S094 Digest Buffer	21
S097 Pre-Wetting Solution	53
S098 Spotting Solution	88
S099 QuantiBlot Wash Solution	77
S100 QuantiBlot DNA Standards	75
S101 SDS, 0.1%	80
S102 Heparin Solution	35
S104 Chromogen Solution	19
S107 TNE, 1X	99
S108 SP, 2X	87
S109 PM PCR Reaction Mixture	49
S110 EDTA, 200 mM	30
S111 QUAD STR/PCR Reaction Mixture	73
S112 dNTP's (10mM stock)	42
S113 Acrylamide Solution, 40%	6
S114 Urea (25 g Aliquot)	101
S114a (18 g Aliquot-377 Sequencer)	101A
S114b Urea Urea (10.8 g Aliquot-377 Sequencer)	101B
S115 Ammonium Persulfate (0.5g Aliquot)	7
S116 GS500 + Loading Buffer (1:3)	33
S117 BSA Solution, 5 mg/mL	9
S118 Primer, VWA 1 (50 $\mu$ M)	54
S119 Primer, VWA 2 (50 $\mu$ M)	56
S120 Primer THO1 1 (50 $\mu$ M)	58
S121 Primer THO1 2 (50 $\mu$ M)	60
S122 Primer F13A1 1 (50 $\mu$ M)	62
S123 Primer, F13A1 2 (50 $\mu$ M)	64
S124 Primer, FES/FPS 1 (50 $\mu$ M)	66
S125 Primer FES/FPS 2 (50 $\mu$ M)	68
S126 Positive Control-External	50
S127 Positive Control-Internal	51
S128 TE <sup>-4</sup> , 1X	94
S131 Blue Formamide + Loading Buffer (5:1)	8A
S132 Sequencing Loading Buffer	79A

Initials: *LC*

Date: *10/21/96*

# Alphabetical Table of Contents

S113 Acrylamide Solution, 40%	6
S115 Ammonium Persulfate (0.5g Aliquot)	7
S018 Analytical Gel Loading Buffer	8
S117 BSA Solution, 5 mg/mL	9
S060 Calibration Control	10
S010 Cell Lysis Buffer (CLB)	12
S064 Cell Pellet Control	13
S022 Chelex, 5%	15
S082 Chelex, 20%	17
S027 Chloroform-Isoamyl Alcohol	18
S104 Chromogen Solution	19
S013 Denaturation Solution	20
S094 Digest Buffer	21
S061 Digestion Control	22
S003 DQ $\alpha$ Citrate Buffer	24
S004 DQ $\alpha$ Hybridization Solution	25
S005 DQ $\alpha$ Wash Solution	26
S093 DTT, 1M	27
S030 DTT, 0.39M	28
S009 EDTA, 0.5M	29
S110 EDTA, 200 mM	30
S056 Ethanol, 70%	31
S031 Gel Neutralization Buffer	32
S116 GS500 + Loading Buffer (1:3)	33
S016 Hae III Buffer, 10X	34
S102 Heparin Solution	35
S026 Hybridization Solution, RFLP	36
S079 Hydrogen Peroxide, 3%	37
S032 Lambda Marker	38
S015 Lithium Chloride, 7.5 M	40
S055 Magnesium Chloride, 0.1M	41
S112 dNTP's (10mM stock)	42
S006 Phenol	43
S053 Phenol-Chloroform-Isoamyl Alcohol	45
S042 Phi-X Marker	46
S034 Phosphate Buffered Saline (PBS)	48
S109 PM PCR Reaction Mixture	49
S126 Positive Control-External	50
S127 Positive Control-Internal	51
S097 Pre-Wetting Solution	53
S118 Primer, VWA 1 (50 $\mu$ M)	54
S119 Primer, VWA 2 (50 $\mu$ M)	56
S120 Primer THO1 1 (50 $\mu$ M)	58
S121 Primer THO1 2 (50 $\mu$ M)	60
S122 Primer F13A1 1 (50 $\mu$ M)	62
S123 Primer, F13A1 2 (50 $\mu$ M)	64

Initials: *RSJ*

Date: *10/7/96*

# Alphabetical Table of Contents

S113 Acrylamide Solution, 40%	6
S115 Ammonium Persulfate (0.5g Aliquot)	7
S018 Analytical Gel Loading Buffer	8
S117 BSA Solution, 5 mg/mL	9
S060 Calibration Control	10
S010 Cell Lysis Buffer (CLB)	12
S064 Cell Pellet Control	13
S022 Chelex, 5%	15
S022 Chelex, 5%	16
S082 Chelex, 20%	17
S027 Chloroform-Isoamyl Alcohol	18
S104 Chromogen Solution	19
S013 Denaturation Solution	20
S094 Digest Buffer	21
S061 Digestion Control	22
S003 DQ $\alpha$ Citrate Buffer	24
S004 DQ $\alpha$ Hybridization Solution	25
S005 DQ $\alpha$ Wash Solution	26
S093 DTT, 1M	27
S030 DTT, 0.39M	28
S009 EDTA, 0.5M	29
S110 EDTA, 200 mM	30
S056 Ethanol, 70%	31
S031 Gel Neutralization Buffer	32
S116 GS500 + Loading Buffer (1:3)	33
S016 Hae III Buffer, 10X	34
S102 Heparin Solution	35
S026 Hybridization Solution, RFLP	36
S079 Hydrogen Peroxide, 3%	37
S032 Lambda Marker	38
S015 Lithium Chloride, 7.5 M	40
S055 Magnesium Chloride, 0.1M	41
S112 dNTP's (10mM stock)	42
S006 Phenol	43
S053 Phenol-Chloroform-Isoamyl Alcohol	45
S042 Phi-X Marker	46
S034 Phosphate Buffered Saline (PBS)	48
S109 PM PCR Reaction Mixture	49
S126 Positive Control-External	50
S127 Positive Control-Internal	51
S097 Pre-Wetting Solution	53
S118 Primer, VWA 1 (50 $\mu$ M)	54
S119 Primer, VWA 2 (50 $\mu$ M)	56
S120 Primer THO1 1 (50 $\mu$ M)	58
S121 Primer THO1 2 (50 $\mu$ M)	60
S122 Primer F13A1 1 (50 $\mu$ M)	62

Initials: RC

Date: 10/20/96

S124 Primer, FES/FPS 1 (50 $\mu$ M)	66
S125 Primer FES/FPS 2 (50 $\mu$ M)	68
S011 Protein Lysis Buffer (PLB)	70
S014 Proteinase-K Enzyme, 10mg/mL	71
S037 Proteinase-K Enzyme, 20mg/mL	72
S111 QUAD STR/PCR Reaction Mixture	73
S100 QuantiBlot DNA Standards	75
S099 QuantiBlot Wash Solution	77
S054 Sarkosyl, 10%	78
S040 Sarkosyl, 20%	79
S101 SDS, 0.1%	80
S045 SDS, 10%	81
S001 SDS, 20%	82
S046 SLS, 20%	83
S043 Sodium Acetate, 2 M	84
S044 Sodium Acetate, 0.2M	85
S035 SP, 25X	86
S108 SP, 2X	87
S098 Spotting Solution	88
S002 SSPE, 20X	89
S047 Stain Extraction Buffer	90
S059 Sterile Water	91
S023 Stripping Solution	92
S017 TRIS-HCl, 1M - PH 8.0	93
S128 TE <sup>-4</sup> , 1X	94
S050 Test Gel Loading Buffer	95
S062 Test Gel Standard	96
S051 TNE, 10X	98
S107 TNE, 1X	99
S052 TRIS, 0.1M - pH 7.8	100
S114 Urea (25 g Aliquot)	101
S024 Wash Solution #1	102
S025 Wash Solution #2	103
S020 Yield Calibrators	104
S021 Yield Gel Loading Buffer	106



Initials: *RCJ*

Date: *10/7/96*

S123 Primer, F13A1 2 (50 $\mu$ M) . . . . .	64
S124 Primer, FES/FPS 1 (50 $\mu$ M) . . . . .	66
S125 Primer FES/FPS 2 (50 $\mu$ M) . . . . .	68
S011 Protein Lysis Buffer (PLB) . . . . .	70
S014 Proteinase-K Enzyme, 10mg/mL . . . . .	71
S037 Proteinase-K Enzyme, 20mg/mL . . . . .	72
S111 QUAD STR/PCR Reaction Mixture . . . . .	73
S100 QuantiBlot DNA Standards . . . . .	75
S099 QuantiBlot Wash Solution . . . . .	77
S054 Sarkosyl, 10% . . . . .	78
S040 Sarkosyl, 20% . . . . .	79
S101 SDS, 0.1% . . . . .	80
S045 SDS, 10% . . . . .	81
S001 SDS, 20% . . . . .	82
S046 SLS, 20% . . . . .	83
S043 Sodium Acetate, 2 M . . . . .	84
S044 Sodium Acetate, 0.2M . . . . .	85
S035 SP, 25X . . . . .	86
S108 SP, 2X . . . . .	87
S098 Spotting Solution . . . . .	88
S002 SSPE, 20X . . . . .	89
S047 Stain Extraction Buffer . . . . .	90
S059 Sterile Water . . . . .	91
S023 Stripping Solution . . . . .	92
S039 TE, 1X . . . . .	93
S049 TE, 100X . . . . .	94
S050 Test Gel Loading Buffer . . . . .	95
S062 Test Gel Standard . . . . .	96
S051 TNE, 10X . . . . .	98
S107 TNE, 1X . . . . .	99
S052 TRIS, 0.1M - pH 7.8 . . . . .	100
S114 Urea (25 g Aliquot) . . . . .	101
S024 Wash Solution #1 . . . . .	102
S025 Wash Solution #2 . . . . .	103
S020 Yield Calibrators . . . . .	104
S021 Yield Gel Loading Buffer . . . . .	106

Initials: *RCJ*

Date: *10/21/96*

# Numerical Table of Contents

S001 SDS, 20%	82
S002 SSPE, 20X	89
S003 DQα Citrate Buffer	24
S004 DQα Hybridization Solution	25
S005 DQα Wash Solution	26
S006 Phenol	43
S009 EDTA, 0.5M	29
S010 Cell Lysis Buffer (CLB)	12
S011 Protein Lysis Buffer (PLB)	70
S013 Denaturation Solution	20
S014 Proteinase-K Enzyme, 10mg/mL	71
S015 Lithium Chloride, 7.5 M	40
S016 Hae III Buffer, 10X	34
S017 TRIS-HCl, 1M - PH 8.0	95
S018 Analytical Gel Loading Buffer	8
S020 Yield Calibrators	104
S021 Yield Gel Loading Buffer	106
S022 Chelex, 5%	15
S023 Stripping Solution	92
S024 Wash Solution #1	102
S025 Wash Solution #2	103
S026 Hybridization Solution, RFLP	36
S027 Chloroform-Isoamyl Alcohol	18
S030 DTT, 0.39M	28
S031 Gel Neutralization Buffer	32
S032 Lambda Marker	38
S034 Phosphate Buffered Saline (PBS)	48
S035 SP, 25X	86
S037 Proteinase-K Enzyme, 20mg/mL	72
S040 Sarkosyl, 20%	79
S042 Phi-X Marker	46
S043 Sodium Acetate, 2 M	84
S044 Sodium Acetate, 0.2M	85
S045 SDS, 10%	81
S046 SLS, 20%	83
S047 Stain Extraction Buffer	90
S050 Test Gel Loading Buffer	95
S051 TNE, 10X	98
S052 TRIS, 0.1M - pH 7.8	100
S053 Phenol-Chloroform-Isoamyl Alcohol	45
S054 Sarkosyl, 10%	78
S055 Magnesium Chloride, 0.1M	41
S056 Ethanol, 70%	31
S059 Sterile Water	91
S060 Calibration Control	10
S061 Digestion Control	22

Initials: *RSJ*

Date: *10/1/96*

# Numerical Table of Contents

S001 SDS, 20%	80
S002 SSPE, 20X	87
S003 DQ $\alpha$ Citrate Buffer	24
S004 DQ $\alpha$ Hybridization Solution	25
S005 DQ $\alpha$ Wash Solution	26
S006 Phenol	43
S009 EDTA, 0.5M	29
S010 Cell Lysis Buffer (CLB)	12
S011 Protein Lysis Buffer (PLB)	68
S013 Denaturation Solution	20
S014 Proteinase-K Enzyme, 10mg/mL	69
S015 Lithium Chloride, 7.5 M	40
S016 Hae III Buffer, 10X	34
S018 Analytical Gel Loading Buffer	8
S020 Yield Calibrators	102
S021 Yield Gel Loading Buffer	104
S022 Chelex, 5%	16
S022 Chelex, 5%	15
S023 Stripping Solution	90
S024 Wash Solution #1	100
S025 Wash Solution #2	101
S026 Hybridization Solution, RFLP	36
S027 Chloroform-Isoamyl Alcohol	18
S030 DTT, 0.39M	28
S031 Gel Neutralization Buffer	32
S032 Lambda Marker	38
S034 Phosphate Buffered Saline (PBS)	48
S035 SP, 25X	84
S037 Proteinase-K Enzyme, 20mg/mL	70
S039 TE, 1X	91
S040 Sarkosyl, 20%	77
S042 Phi-X Marker	46
S043 Sodium Acetate, 2 M	82
S044 Sodium Acetate, 0.2M	83
S045 SDS, 10%	79
S046 SLS, 20%	81
S047 Stain Extraction Buffer	88
S049 TE, 100X	92
S050 Test Gel Loading Buffer	93
S051 TNE, 10X	96
S052 TRIS, 0.1M - pH 7.8	98
S053 Phenol-Chloroform-Isoamyl Alcohol	45
S054 Sarkosyl, 10%	76
S055 Magnesium Chloride, 0.1M	41
S056 Ethanol, 70%	31
S059 Sterile Water	89

Initials: *RCI* Date: *10/21/96*

S062 Test Gel Standard	96
S064 Cell Pellet Control	13
S079 Hydrogen Peroxide, 3%	37
S082 Chelex, 20%	17
S093 DTT, 1M	27
S094 Digest Buffer	21
S097 Pre-Wetting Solution	53
S098 Spotting Solution	88
S099 QuantiBlot Wash Solution	77
S100 QuantiBlot DNA Standards	75
S101 SDS, 0.1%	80
S102 Heparin Solution	35
S104 Chromogen Solution	19
S107 TNE, 1X	99
S108 SP, 2X	87
S109 PM PCR Reaction Mixture	49
S110 EDTA, 200 mM	30
S111 QUAD STR/PCR Reaction Mixture	73
S112 dNTP's (10mM stock)	42
S113 Acrylamide Solution, 40%	6
S114 Urea (25 g Aliquot)	101
S115 Ammonium Persulfate (0.5g Aliquot)	7
S116 GS500 + Loading Buffer (1:3)	33
S117 BSA Solution, 5 mg/mL	9
S118 Primer, VWA 1 (50 $\mu$ M)	54
S119 Primer, VWA 2 (50 $\mu$ M)	56
S120 Primer THO1 1 (50 $\mu$ M)	58
S121 Primer THO1 2 (50 $\mu$ M)	60
S122 Primer F13A1 1 (50 $\mu$ M)	62
S123 Primer, F13A1 2 (50 $\mu$ M)	64
S124 Primer, FES/FPS 1 (50 $\mu$ M)	66
S125 Primer FES/FPS 2 (50 $\mu$ M)	68
S126 Positive Control-External	50
S127 Positive Control-Internal	51
S128 TE <sup>-4</sup> , 1X	94

Initials: *RS*

Date: *10/7/96*

S060 Calibration Control	10
S061 Digestion Control	22
S062 Test Gel Standard	94
S064 Cell Pellet Control	13
S079 Hydrogen Peroxide, 3%	37
S082 Chelex, 20%	17
S093 DTT, 1M	27
S094 Digest Buffer	21
S097 Pre-Wetting Solution	51
S098 Spotting Solution	86
S099 QuantiBlot Wash Solution	75
S100 QuantiBlot DNA Standards	73
S101 SDS, 0.1%	78
S102 Heparin Solution	35
S104 Chromogen Solution	19
S107 TNE, 1X	97
S108 SP, 2X	85
S109 PM PCR Reaction Mixture	49
S110 EDTA, 200 mM	30
S111 QUAD STR/PCR Reaction Mixture	71
S112 dNTP's (10mM stock)	42
S113 Acrylamide Solution, 40%	6
S114 Urea (25 g Aliquot)	99
S115 Ammonium Persulfate (0.5g Aliquot)	7
S116 GS500 + Loading Buffer (1:3)	33
S117 BSA Solution, 5 mg/mL	9
S118 Primer, VWA 1 (50 $\mu$ M)	52
S119 Primer, VWA 2 (50 $\mu$ M)	54
S120 Primer THO1 1 (50 $\mu$ M)	56
S121 Primer THO1 2 (50 $\mu$ M)	58
S122 Primer F13A1 1 (50 $\mu$ M)	60
S123 Primer, F13A1 2 (50 $\mu$ M)	62
S124 Primer, FES/FPS 1 (50 $\mu$ M)	64
S125 Primer FES/FPS 2 (50 $\mu$ M)	66
S126 Positive Control-External	50
S127 Positive Control-Internal	51

Initials: RSDate: 10/7/96**S113 Acrylamide Solution, 40%**

lot number: \_\_\_\_\_

standard batch sizes: 250 mL, 375 mL

**Ingredients**

	amount	amount
RM2085 polyacrylamide bisacrylamide (19:1) (premixed)	100 $\pm$ 1 g	150 $\pm$ 1.5 g
RM0409 mixed bed ion exchange resin	20 $\pm$ 1 g	30 $\pm$ 1.5 g
<u>Initial amount deionized water</u>	<u>150 <math>\pm</math> 5 mL</u>	<u>225 <math>\pm</math> 5 mL</u>
Final Volume	250 $\pm$ 5 mL	375 $\pm$ 5 mL

**Procedure**

Add initial amount deionized water to the original polyacrylamide bisacrylamide bottle.

Stir at low heat (position 2 on hot plate) until dissolved.

Adjust to the final volume with deionized water.

Pour back into bottle.

Add mixed bed ion exchange resin and stir for 30 minutes.

**Remove the stir bar.**

Store at 2-8°C for one month. Label with expiration date.

**Data Log**

	source	lot	amount
RM2085 Polyacrylamide Bisacrylamide	_____	_____	_____
RM0409 mixed bed ion exchange resin	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RG

Date: 10/7/96

**S115 Ammonium Persulfate (0.5g Aliquot)**

lot number: \_\_\_\_\_

standard batch size: ~ 25 tubes x 0.5g

Ingredients	aliquot	total amount
RM 0016 ammonium persulfate electrophoresis grade	$0.5 \pm 0.05$ g	$12.5 \pm 1$ g

**Procedure**

**NOTE: WHEN WORKING WITH POWDERED AMMONIUM PERSULFATE WEAR GLOVES, EYE PROTECTION, LAB COAT, AND RESPIRATOR FOR SAFETY.**

Fill out chemical logbook.

Using weigh paper, weigh  $0.5 \pm 0.05$  g aliquots of ammonium persulfate.

Transfer the aliquots to 15mL conical tubes.

Cap all tubes tightly and label rack containing tubes with contents, lot number, date, initials, and safety data.

Store at room temperature.

**Data Log**

source                      lot                      amount

RM0016 ammonium persulfate      \_\_\_\_\_      \_\_\_\_\_      \_\_\_\_\_

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCIDate: 10/2/96**S018 Analytical Gel Loading Buffer** lot number: \_\_\_\_\_

standard batch size: 100 mL

Ingredients	final concentration	amount
RM020 bromophenol blue	0.25%	0.25 ± 0.01 g
RM217 xylene cyanol	0.25%	0.25 ± 0.01 g
RM040 Ficoll 400	12.5%	12.5 ± 0.1 g
S009 EDTA, 0.5M	50. mM	10.0 ± 0.1 mL
RM083 TAE, 10X	5.0 X	50.0 ± 0.5 mL

**Procedure**

Combine the TAE, EDTA, and Ficoll.

Mix well. The solution may need to be heated gently to dissolve the Ficoll.

Add the bromophenol blue and xylene cyanol.

Mix well.

When all the solids are dissolved, bring up to volume using deionized water.

Filter sterilize.

Dispense 1.5 mL aliquots into 1.5 mL eppendorf tubes.

Store at -20°C.

Data Log	source	lot	amount
RM020 bromophenol blue	_____	_____	_____
RM217 xylene cyanol	_____	_____	_____
RM040 Ficoll 400	_____	_____	_____
S009 EDTA, 0.5M	_____	_____	_____
RM083 TAE, 10X	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_



Initials: RCJ

Date: 2/9/98

**S131 Blue Formamide + Loading Buffer (5:1)**

lot number: \_\_\_\_\_

standard batch size: 20 x 1200  $\mu$ L

Ingredients	amount	final ratio
RM2099 Formamide	1000 $\pm$ 20 $\mu$ L	5
S132 Sequencing Loading Buffer	200 $\pm$ 10 $\mu$ L	1

**Procedure**

**NOTE: PREPARE AWAY FROM AMPLIFIED DNA TO MINIMIZE CONTAMINATION. USING CLEAN GLOVES IS ESSENTIAL; CHANGE THEM AS OFTEN AS NEEDED.**

Clean the bench top thoroughly using a 10% bleach solution, and cover it with new bench paper.

Label 20 1.5mL reaction tubes.

Add Formamide to each tube. Add blue sequencing buffer to each tube.

Close all tubes and mix.

Store at 2-8°C.

**Data Log**

	source	lot	amount
RM2099 Formamide	_____	_____	_____
S132 Sequencing Loading Buffer	_____	_____	_____

**Quality Control**

QC032 STR gel electrophoresis

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RC

Date: 10/7/21

**S117 BSA Solution, 5 mg/mL**

lot number: \_\_\_\_\_

standard batch size: 25 mL

**Ingredients**

amount

RM222 Bovine serum albumin

125 mg

S059 Sterile water

25 mL (guideline)

**Procedure**

Autoclave a 50 mL glass beaker with a stir bar in it.

Add the BSA to 20 mL of sterile water in the glass beaker.

Stir gently over **very low** heat until the BSA is completely dissolved.

Add the solution to a 50 mL disposable conical tube.

Add sterile water to a final volume of 25 mL.

Aliquot approximately 0.5 mL of BSA solution into 1.5 mL microcentrifuge tubes.

Label each tube with "BSA" and the lot number. Label the rack with expiration date.

Store at -20°C.

**Data Log**

source

lot

amount

RM222 BSA

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

S059 Sterile water

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**Quality Control**

QC023 QuantiBlot Quality Control of Solutions- test 20 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *RCJ*Date: *10/7/96***S060 Calibration Control**

lot number: \_\_\_\_\_

page 1 of 2

Ingredients	initial concentration (ng/μL)	initial volume (μL)	final concentration	final volume (μL)
RM221 K562 DNA			7.5 ng/μL	
S021 yield gel loading buffer	5 X		1 X	----
S059 sterile water	-----		----	----

**Calculations**

Record the initial concentration in ng/μL and the initial volume in μL of the K562 DNA received from the manufacturer.

Calculate the final volume according to equation 1.

$$(\text{final volume}) = \frac{(\text{initial DNA concentration})(\text{initial DNA volume})}{(7.5 \text{ ng/}\mu\text{L})} \quad \text{equation 1}$$

Record the final volume above. The final volume is the total batch size.

Calculate the amount of buffer to be added according to equation 2.

$$(\text{buffer volume}) = 0.2(\text{final volume}) \quad \text{equation 2}$$

Calculate the amount of sterile water to be added according to equation 3.

$$(\text{water volume}) = [0.8 * (\text{final volume})] - (\text{initial DNA volume}) \quad \text{equation 3}$$

Record the buffer and water volumes above.

To check the calculations, add together the initial volumes of DNA, loading buffer, and sterile water.

The sum of the initial volumes must be equal to the calculated final volume.

Initials: RCJ

Date: 10/7/96

**S060 Calibration Control**

lot number: \_\_\_\_\_

page 2 of 2

**Procedure**

Combine the DNA, loading buffer, and sterile water.

Mix well.

Using sterile pipet tips, dispense 200  $\mu$ L aliquots into sterile 1.5 mL eppendorf tubes.

Store at -20°C.

**Data Log**

source	lot	amount
--------	-----	--------

RM221 K562 DNA	_____	_____
----------------	-------	-------

S021 yield gel loading buffer	_____	_____
-------------------------------	-------	-------

S059 sterile water	_____	_____
--------------------	-------	-------

**Quality Control**

QC026 Gel Electrophoresis

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCJDate: 10/7/96**S010 Cell Lysis Buffer (CLB)**

lot number: \_\_\_\_\_

standard batch size: 2 L

Ingredients	final concentration	amount
RM068 sucrose	320 mM	219 ± 3 g
RM073 TRIS	10 mM	2.4 ± 0.1 g
RM046 magnesium chloride, hexahydrate	5 mM	2.0 ± 0.1 g
RM075 Triton X-100	1.0 %	20 ± 1 mL
RM096 hydrochloric acid	---	

**Procedure**

Dissolve the sucrose, TRIS, and magnesium chloride in approximately 1.5 L deionized water.

Add the Triton to the solution.

Adjust the pH to 7.6 with hydrochloric acid

Mix well.

Adjust the volume to 2 L with deionized water.

Filter sterilize.

Dispense into sterile 50 mL centrifuge tubes.

Store at 2-8°C.

Data Log	source	lot	amount
RM068 sucrose	_____	_____	_____
RM073 TRIS	_____	_____	_____
RM046 magnesium chloride, hexahydrate	_____	_____	_____
RM075 Triton X-100	_____	_____	_____
RM096 hydrochloric acid	_____	_____	_____

**Quality Control**

QC023 QuantiBlot Quality Control of Solutions- test 150 µL of solution

final pH: \_\_\_\_\_ spec: 7.6 ± 0.1

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *RG*Date: *10/17/96***S064 Cell Pellet Control**

lot number: \_\_\_\_\_

page 1 of 2

Ingredients	concentration of cells	total volume (mL)	cells per aliquot	aliquot volume (mL)
RM243 K562 cells			$1 \times 10^6$	
S034 phosphate buffered saline (PBS)	----	----	----	----

**Calculations**

Record the concentration of K562 cells in the suspension received from the manufacturer.

Record the total volume. This is the batch size.

Calculate the volume (in mL) which yields  $1 \times 10^6$  cells according to equation 1.

$$(\text{aliquot volume}) = \frac{(1 \times 10^6 \text{ cells})}{(\text{concentration of cells})} \quad \text{equation 1}$$

The aliquot volume must fit into a 1.5 mL eppendorf tube. The concentration of the cell suspension may have to be adjusted.

If the cell concentration is too low, the cells may be spun at  $180 \times g$  for 5 minutes at  $4^\circ\text{C}$ . Remove the excess media to give the desired concentration.

If the cell concentration is too high, PBS may be added to reach the desired concentration. After adding PBS, make sure the cells are well suspended before aliquoting.

Record the calculated aliquot volume.

**Procedure**

The following steps must be done on ice or at  $4^\circ\text{C}$ .

Bring the cell suspension up to the desired final volume.

Suspend the cells evenly by pipetting up and down or by gently inverting the container.

Initials: 25

Date: 10/7/88

**S064 Cell Pellet Control**

lot number: \_\_\_\_\_

page 2 of 2

Add aliquots of cell suspension to 1.5 mL eppendorf tubes.

Spin the tubes at 180 x g for 1 minute at 4°C, and remove the excess supernatant.

The tubes can be aliquoted and spun in sets of 52. Each set should be packaged separately in a seal-a-meal bag, labeled with the lot number and numbered sequentially.

Store the bags at -70°C.

**Data Log**

source	lot	amount
--------	-----	--------

RM243 K562 cells	_____	_____
------------------	-------	-------

S034 phosphate buffered saline	_____	_____
--------------------------------	-------	-------

**Quality Control**

QC024 Non-Organic Extraction

QC027 Southern Blotting and Hybridization

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RS

Date: 10/7/96

**S022 Chelex, 5%**

lot number: \_\_\_\_\_

standard batch size: 500 mL

**Ingredients**

final  
concentration

amount

RM027 Chelex 100

5. %

25 ± 2 g

S059 sterile water

---

450 ± 50 mL (guideline)

**Procedure**

Filter sterilize approximately 600 mL deionized water.

Pour the water into a 500 mL bottle.

Save the bottom container from the disposable filter unit.

Autoclave the water at 250°F for 30 minutes.

Add the Chelex to the bottom container of the filter unit.

Allow the water to cool after autoclaving.

Add sterile water to the Chelex to a volume of 500 mL using the graduation markings on the disposable filter container.

Mix on a magnetic stir plate.

While the stock solution is mixing, aliquot 10 mL each into 15 mL centrifuge tubes.

Store at 2-8°C.

**Data Log**

source

lot

amount

RM027 Chelex 100

\_\_\_\_\_

S059 sterile water

\_\_\_\_\_

**Quality Control**

QC014 Chelex Extraction

made by: \_\_\_\_\_

date: \_\_\_\_\_



Initials: RS

Date: 10/7/96

**S022 Chelex, 5%**

lot number: \_\_\_\_\_

standard batch size: 400 mL

Ingredients	final concentration	amount
RM536 DNA Extraction Reagent	5 %	100 mL
S059 sterile water	---	300 mL

**Procedure**

Filter sterilize approximately 300 mL deionized water.

Pour the water into a 500 mL bottle.

Save the bottom container from the disposable filter unit.

Autoclave the water at 250°F for 30 minutes.

Add 100 mL of DNA Extraction Reagent to the bottom container of the filter unit.

Allow the water to cool after autoclaving.

Add sterile water to the the DNA Extraction Reagent to a volume of 400 mL using the graduation markings on the disposable filter container.

Mix on a magnetic stir plate.

While the stock solution is mixing, aliquot 10 mL each into 15 mL centrifuge tubes.

Store at 2-8°C.

Data Log	source	lot	amount
RM536 DNA Extraction Reagent	_____	_____	_____
S059 sterile water	_____	_____	_____

**Quality Control**

QC014 Chelex Extraction

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials:   25  

Date:   10/7/96  

**S082 Chelex, 20%**

lot number: \_\_\_\_\_

standard batch size: 500 mL

Ingredients	final concentration	amount
RM027 Chelex 100	20. %	100 ± 2 g
S059 sterile water	---	450 ± 50 mL (guideline)

### Procedure

Filter sterilize approximately 600 mL deionized water.

Pour the water into a 500 mL bottle.

Save the bottom container from the disposable filter unit.

Autoclave the water at 250°F for 30 minutes.

Add the Chelex to the bottom container of the filter unit.

Allow the water to cool after autoclaving.

Add sterile water to the Chelex to a volume of 500 mL using the graduation markings on the disposable filter container.

Mix on a magnetic stir plate.

While the stock solution is mixing, aliquot 10 mL each into 15 mL centrifuge tubes.

Store at 2-8°C.

Data Log	source	lot	amount
RM027 Chelex 100	_____	_____	_____
S059 sterile water	_____	_____	_____

### Quality Control

QC017 Differential Extraction

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RU

Date: 10/7/96

**S027 Chloroform-Isoamyl Alcohol**

lot number: \_\_\_\_\_

standard batch size: 500 mL

Ingredients	final concentration	amount
RM088 chloroform	96. %	480 $\pm$ 3 mL
RM089 isoamyl alcohol	4. %	20 $\pm$ 3 mL

**Procedure**

**NOTE:** Use only glass graduated cylinders and containers.

Measure the isoamyl alcohol into a 500 mL brown bottle.

Add the chloroform.

Store at 2-8°C in a flammable materials refrigerator.

Data Log	source	lot	amount
RM088 chloroform	_____	_____	_____
RM089 isoamyl alcohol	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials:   AG  

Date:   10/7/96  

**S104 Chromogen Solution**

lot number: \_\_\_\_\_

standard batch size: 30 mL

Ingredients	final concentration	amount
RM435 chromogen:TMB	_____	60 mg
RM009 ethanol, 100% reagent grade	_____	30 mL

**Procedure**

Bring bottle of chromogen:TMB to room temperature.

Before opening, lightly tap the bottle on the counter to bring its contents to the bottom.

Carefully remove the stopper and reconstitute the chromogen:TMB with the room temperature ethanol.

**CAUTION: DO NOT USE ETHANOL STORED IN A METAL CONTAINER;  
ONLY USE 100% REAGENT GRADE ETHANOL.**

Recap the bottle and seal with Parafilm.

Tilt the bottle several times to ensure that all the powder is removed from within the rubber cap.

Shake on an orbital shaker for about 30 minutes.

Store at 2-8°C and away from rust.

The solution is stable for six months.

Data Log	source	lot	amount
RM435 chromogen	_____	_____	_____
RM009 ethanol, 100%	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: 129

Date: 10/2/96

**S013 Denaturation Solution**

lot number: \_\_\_\_\_

standard batch size: 20 L

**Ingredients**

final  
concentration

amount

RM005 sodium chloride

1.5 M

1750 ± 50 g

RM004 sodium hydroxide, 10N

0.5 M

1000 ± 100 mL

**Procedure**

Measure the sodium chloride into a 20 L carboy with approximately 5-10 L deionized water.

Add the sodium hydroxide solution.

Mix well on magnetic stir plate using a stir bar.

Raise to the final volume with deionized water.

Mix well.

Measure and record the pH using pH paper.

Store at room temperature.

**Data Log**

source

lot

amount

RM005 sodium chloride

\_\_\_\_\_

RM004 sodium hydroxide, 10N

\_\_\_\_\_

**Quality Control**

final pH: \_\_\_\_\_ spec: ≥ 12

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials:   ZC  Date:   10/7/96  **S094 Digest Buffer**

lot number: \_\_\_\_\_

standard batch size: 2 L

Ingredients	final concentration	amount
S009 EDTA, 0.5M	10. mM	40 ± 2 mL
RM073 TRIS	10. mM	2.4 ± 0.2 g
RM0005 sodium chloride	50. mM	5.8 ± 0.4 g
S001 SDS, 20%	2.0 %	200 ± 2 mL
RM096 hydrochloric acid	---	---

**Procedure**

Add the EDTA, TRIS, sodium chloride, and SDS to approximately 1.5 L deionized water.  
Adjust the pH to 7.5 with hydrochloric acid.

Bring up to the final volume with deionized water.

Mix well.

Measure and record the final pH.

Aliquot into 50 mL centrifuge tubes.

Store at room temperature.

Data Log	source	lot	amount
S009 EDTA, 0.5M	_____	_____	_____
RM073 TRIS	_____	_____	_____
RM0005 sodium chloride	_____	_____	_____
S001 SDS, 20%	_____	_____	_____
RM096 hydrochloric acid	_____	_____	_____

**Quality Control**

final pH: \_\_\_\_\_ specification: 7.5 ± 0.1

QC023 QuantiBlot Quality Control of Solutions- Test 150 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *Re*

Date: *10/4/96*

**S061 Digestion Control**

lot number: \_\_\_\_\_

page 1 of 2

INGREDIENTS	initial concentration (ng/ $\mu$ L)	initial volume ( $\mu$ L)	final concentration	final volume ( $\mu$ L)
RM221 K562 DNA			2 ng/ $\mu$ L	
S128 TE <sup>-4</sup> , 1X	1 X		----	----

**Calculations**

Record the initial concentration in ng/ $\mu$ L and the initial volume in  $\mu$ L of the K562 DNA received from the manufacturer.

Calculate the final volume according to equation 1.

$$(\text{final volume}) = \frac{(\text{initial DNA concentration})(\text{initial DNA volume})}{(2 \text{ ng}/\mu\text{L})} \quad \text{equation 1}$$

Record the final volume above. The final volume is the total batch size.

Calculate the amount of 1X TE<sup>-4</sup> to be added according to equation 2.

$$(\text{TE}^{-4} \text{ volume}) = (\text{final volume}) - (\text{initial DNA volume}) \quad \text{equation 2}$$

Record the TE<sup>-4</sup> volume above.

To check the calculations, add together the initial volumes of DNA and 1X TE<sup>-4</sup>.

The sum of the initial volumes must be equal to the calculated final volume.

Initials:

Date:

**S061 Digestion Control**

lot number: \_\_\_\_\_

page 2 of 2

**Procedure**

Combine the DNA and 1X TE<sup>-4</sup>. **Make sure the DNA is in solution. If the DNA is frozen, resuspend for at least 2 hr RT or 4°C overnight**

Mix well.

Using sterile pipet tips, dispense 250 µL aliquots into 1.8 mL microcentrifuge tubes.

Store at -20°C.

**Data Log**

source

lot

amount

RM221 K562 DNA

\_\_\_\_\_

S128 TE<sup>-4</sup>, 1X

\_\_\_\_\_

**Quality Control**

QC026 Gel Electrophoresis

made by: \_\_\_\_\_ date: \_\_\_\_\_



Initials: RL

Date: 10/7/86

**S061 Digestion Control**

lot number: \_\_\_\_\_

page 1 of 2

INGREDIENTS	initial concentration (ng/μL)	initial volume (μL)	final concentration	final volume (μL)
RM221 K562 DNA			2 ng/μL	
S039 TE, 1X	1 X		----	----

**Calculations**

Record the initial concentration in ng/μL and the initial volume in μL of the K562 DNA received from the manufacturer.

Calculate the final volume according to equation 1.

$$(\text{final volume}) = \frac{(\text{initial DNA concentration})(\text{initial DNA volume})}{(2 \text{ ng/}\mu\text{L})} \quad \text{equation 1}$$

Record the final volume above. The final volume is the total batch size.

Calculate the amount of 1X TE to be added according to equation 2.

$$(\text{TE volume}) = (\text{final volume}) - (\text{initial DNA volume}) \quad \text{equation 2}$$

Record the TE volume above.

To check the calculations, add together the initial volumes of DNA and 1X TE.

The sum of the initial volumes must be equal to the calculated final volume.

Initials: RL

Date: 10/7/96

**S061 Digestion Control**

lot number: \_\_\_\_\_

page 2 of 2

**Procedure**

Combine the DNA and 1X TE. **Make sure the DNA is in solution. If the DNA is frozen, resuspend for at least 2 hr RT or 4°C overnight**

Mix well.

Using sterile pipet tips, dispense 250 µL aliquots into 1.8 mL microcentrifuge tubes.

Store at -20°C.

**Data Log**

source	lot	amount
--------	-----	--------

RM221 K562 DNA	_____	_____
----------------	-------	-------

S039 TE, 1X	_____	_____
-------------	-------	-------

**Quality Control**

QC026 Gel Electrophoresis

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RS

Date: 10/7/96

**S003 DQ $\alpha$  Citrate Buffer**

lot number: \_\_\_\_\_

standard batch size: 4 L

**Ingredients**

final  
concentration

amount

RM001 trisodium citrate

-----

73.6  $\pm$  0.1 g

RM002 citric acid

-----

24  $\pm$  1 g (guideline)

**Procedure**

Dissolve the sodium citrate in approximately 3 liters deionized water.

Adjust the pH to 5.0 by addition of citric acid (approximately 24 g).

Adjust the final volume to 4 liters with deionized water.

Mix well.

Measure and record the final pH.

Dispense into a 4 L bottle.

Store at room temperature.

**Data Log**

source

lot

amount

RM001 trisodium citrate

\_\_\_\_\_

RM002 citric acid

\_\_\_\_\_

**Quality Control**

final pH: \_\_\_\_\_ specification 5.0  $\pm$  0.2

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RSJ

Date: 10/7/96

**S004 DQ $\alpha$  Hybridization Solution**

lot number: \_\_\_\_\_

standard batch size: 4 L

**Ingredients**

final  
concentration

amount

S002 SSPE, 20X

5.0 X

1000  $\pm$  10 mL

S001 SDS, 20%

0.50 %

100  $\pm$  1 mL

**Procedure**

Combine the SSPE and 2.9 L deionized water in a 4 L flask.

Add the SDS.

Warm the solution until all solids are dissolved.

Mix well.

Dispense into 1 L bottles.

Store at room temperature.

**Data Log**

source

lot

amount

S002 SSPE, 20X

\_\_\_\_\_

S001 SDS, 20%

\_\_\_\_\_

**Quality Control**

QC016 DQ $\alpha$  Hybridization

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials:   ZC  

Date:   10/2/96  

**S005 DQ $\alpha$  Wash Solution**

lot number: \_\_\_\_\_

standard batch size: 4 L

**Ingredients**

final  
concentration

amount

S002 SSPE, 20X

2.5 X

500  $\pm$  10 mL

S001 SDS, 20%

0.10 %

20  $\pm$  1 mL

**Procedure**

Measure 20 mL 20% SDS in a 50 mL graduated cylinder.

Raise the volume of the SDS solution to 50 mL by adding 30 mL deionized water.

Pour the SDS into a 4 L bottle.

Add 500 mL SSPE and 3450 mL deionized water.

Cap and mix well by inverting.

Store at room temperature.

**Data Log**

source

lot

amount

S002 SSPE, 20X

\_\_\_\_\_

S001 SDS, 20%

\_\_\_\_\_

**Quality Control**

QC016 DQ $\alpha$  hybridization

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RS

Date: 10/3/96

**S093 DTT, 1M**

lot number: \_\_\_\_\_

standard batch size: 5 mL

**Ingredients**

final  
concentration

amount

RM101 dithiothreitol

1.0 M

0.77 ± 0.05 g

S059 sterile water

-----

-----

**Procedure**

Add the DTT to approximately 4 mL sterile, deionized water in a 15 mL centrifuge tube.

Mix well.

When the DTT is dissolved, bring up to volume with sterile, deionized water.

Filter sterilize.

Dispense 250 µL aliquots into sterile 0.5 mL eppendorf tubes.

Store at -20°C.

**Data Log**

source

lot

amount

RM101 dithiothreitol

\_\_\_\_\_

S059 sterile water

\_\_\_\_\_

**Quality Control**

QC023 QuantiBlot Quality Control of Solutions- Test 20 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: AC

Date: 10/7/96

**S030 DTT, 0.39M**

lot number: \_\_\_\_\_

standard batch size: 1.5 mL

**Ingredients**

final  
concentration

amount

RM101 dithiothreitol

0.39 M

0.090 ± 0.001 g

S059 sterile water

—

1.5 mL (Guideline)

**Procedure**

Add the DTT to approximately 1 mL sterile water in a sterile 1.5 mL eppendorf tube.

Mix well.

When the DTT is dissolved, bring up to volume with sterile water.

Dispense 500 µL aliquots into sterile 0.5 mL eppendorf tubes.

Store at -20°C.

Discard after 6 months

**Data Log**

source

lot

amount

RM101 dithiothreitol

\_\_\_\_\_

S059 sterile water

\_\_\_\_\_

**Quality Control**

QC023 QuantiBlot Quality Control of Solutions- Test 40 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RSJ

Date: 10/2/96

**S009 EDTA, 0.5M**

lot number: \_\_\_\_\_

standard batch size: 500 mL

Ingredients	final concentration	amount
RM003 EDTA	0.50 M	93 ± 1 g
RM004 sodium hydroxide, 10N	-----	-----

**Procedure**

Add the EDTA to approximately 250 mL deionized water.

Adjust the pH to 8.0 with sodium hydroxide solution.

Mix well.

When the EDTA is dissolved, adjust the pH to 8.0.

Bring up to volume with deionized water.

Check and record the final pH.

Dispense into 125 mL bottles.

Autoclave at 250°F for 20 minutes.

Store at room temperature.

**Data Log**

	source	lot	amount
RM003 EDTA	_____	_____	_____
RM004 sodium hydroxide, 10N	_____	_____	_____

**Quality Control**

final pH: \_\_\_\_\_ specification: 8.0 ± 0.1

made by: \_\_\_\_\_ date: \_\_\_\_\_



Initials: RU

Date: 10/19/96

**S110 EDTA, 200 mM**

lot number: \_\_\_\_\_

standard batch size: 25 mL

**Ingredients**

final  
concentration

amount

S009 EDTA, 0.5 M

200 mM

10.0 ± 0.5 mL

**Procedure**

Add the EDTA to 15 mL deionized water.

Mix thoroughly.

Filter sterilize or autoclave at 250°F for 20 minutes.

Dispense into 300 µL aliquots.

Store at room temperature.

**Data Log**

source

lot

amount

S009 EDTA, 0.5 M

\_\_\_\_\_

**Quality Control**

QC016 PCR Kit Hybridization- Test 20 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RU

Date: 10/7/96

**S056 Ethanol, 70%**

lot number: \_\_\_\_\_

standard batch size: 500 mL

**Ingredients**

final  
concentration

amount

RM009 ethanol, 100%

70%

350 ± 10 mL

**Procedure**

Measure the ethanol in a graduated cylinder.

Bring up to volume with deionized water.

Dispense into a storage container.

Store at room temperature or at 2-8°C.

**Data Log**

source

lot

amount

RM009 ethanol, 100%

\_\_\_\_\_

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RLS

Date: 10/7/96

**S031 Gel Neutralization Buffer**

lot number: \_\_\_\_\_

standard batch size: 8 L

Ingredients	final concentration	amount
RM073 TRIS	0.20 M	194 ± 4 g
RM005 sodium chloride	0.30 M	140 ± 4 g
RM096 hydrochloric acid	-----	100 ± 10 mL (guideline)

**Procedure**

Add the TRIS and NaCl to approximately 6 L deionized water in an 8 L carboy.

Adjust the pH to 7.6 using concentrated HCl.

Bring up to the final volume with deionized water.

Mix well.

Check and record the final pH.

Store at room temperature.

Data Log	source	lot	amount
RM073 TRIS	_____	_____	_____
RM005 sodium chloride	_____	_____	_____
RM096 hydrochloric acid	_____	_____	_____

final pH: \_\_\_\_\_ specification: 7.6 ± 0.1

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCJ

Date: 2/9/98

**S116 GS500 + Loading Buffer (1:3)**

lot number: \_\_\_\_\_

standard batch size: 10 x 800  $\mu$ L

Ingredients	amount	ratio
RM1062 GS 500 Rox labelled	200 $\pm$ 2 $\mu$ L	1
RM0451 Loading Buffer	600 $\pm$ 6 $\mu$ L	3

**Procedure**

**NOTE: PREPARE AWAY FROM AMPLIFIED DNA TO MINIMIZE CONTAMINATION. USING CLEAN GLOVES IS ESSENTIAL; CHANGE THEM AS OFTEN AS NEEDED.**

Clean the bench top thoroughly using a 10% bleach solution, and cover it with new bench paper.

Label a screw cap tube with GS500, add 200 $\mu$ L GS 500 and 600 $\mu$ L loading buffer to the tube.

Close and mix the tube.

Store at 2-8°C.

**Data Log**

	source	lot	amount
RM1062 ABI GS500 Rox labeled	_____	_____	_____
RM0451 ABI Loading Buffer	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: ACJ

Date: 10/17/96

**S116 GS500 + Loading Buffer (1:3)**

lot number: \_\_\_\_\_

standard batch size: 800  $\mu$ L

**Ingredients**

amount

RM1062 GS 500 Rox labelled

200  $\pm$  2  $\mu$ L

RM0451 Loading Buffer

600  $\pm$  6  $\mu$ L

**Procedure**

**NOTE: PREPARE AWAY FROM AMPLIFIED DNA TO MINIMIZE CONTAMINATION. USING CLEAN GLOVES IS ESSENTIAL; CHANGE THEM AS OFTEN AS NEEDED.**

Clean the bench top thoroughly using a 10% bleach solution, and cover it with new bench paper.

Add the loading buffer to the original GS500 tube.

Mark the tube.

Store at 2-8°C.

**Data Log**

	source	lot	amount
RM1062 ABI GS500 Rox labeled	_____	_____	_____
RM0451 ABI Loading Buffer	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: PCDate: 10/7/96**S016 Hae III Buffer, 10X**

lot number: \_\_\_\_\_

standard batch size: 100 mL

Ingredients	final concentration	amount
RM073 TRIS	500 mM	6.0 ± 0.2 g
RM046 magnesium chloride, hexahydrate	100 mM	2.0 ± 0.1 g
RM005 sodium chloride	500 mM	2.9 ± 0.1 g
RM096 hydrochloric acid	---	

**Procedure**

Add the TRIS, magnesium chloride, and sodium chloride to approximately 75 mL deionized water.

Mix well.

Adjust the pH to 8.0 with hydrochloric acid.

Bring up to the final volume with deionized water.

Dispense into a sterile 125 mL bottle.

Autoclave at 250°F for 20 minutes.

Using sterile pipet tips, dispense 1 mL aliquots into sterile 1.5 mL eppendorf tubes.

Store at -20°C.

Data Log	source	lot	amount
RM073 TRIS	_____	_____	_____
RM046 magnesium chloride, hexahydrate	_____	_____	_____
RM005 sodium chloride	_____	_____	_____
RM096 hydrochloric acid	_____	_____	_____

**Quality Control**

final pH: \_\_\_\_\_ spec: 8.0 ± 0.1

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RS

Date: 10/7/86

**S102 Heparin Solution**

lot number: \_\_\_\_\_

standard batch size: 50 mL

Ingredients	final concentration	amount
S035 SP, 25X	5 X	10 ± 0.5 mL
RM028 Heparin	50 mg/mL	2.5 ± 0.05 g
RM061 Na Azide	0.2%	0.1 ± 0.005 g
SO59 Sterile Water	_____	

**Procedure**

Weigh heparin in a sterile 50 mL tube.

Add the Na Azide, SP and mix.

Add the sterile water to volume on tube.

Mix.

The solution may be heated to help dissolve the heparin.

Filter sterilize.

Dispense into a new sterile 50 mL tube.

Store at 4°C.

Data Log	source	lot	amount
S035 SP, 25X	_____	_____	_____
RM028 Heparin	_____	_____	_____
RM061 Na Azide	_____	_____	_____
SO59 Sterile Water	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RC

Date: 10/7/96

**S026 Hybridization Solution, RFLP**

lot number: \_\_\_\_\_

standard batch size: 1 L

Ingredients	final concentration	amount
S035 SP, 25X	1.5 X	60 ± 2 mL
S001 SDS, 20%	7 %	350 ± 10 mL

**Procedure**

Add the SP to 590 mL deionized water.

Add the SDS to the solution.

Heat the solution to dissolve the SDS.

Mix well.

Rinse the filter of a disposable filter unit with approximately 500 mL sterile, deionized water.

Filter sterilize the warm hybridization solution.

Dispense into 250 mL aliquots.

Store at room temperature.

Data Log	source	lot	amount
S035 SP, 25X	_____	_____	_____
S001 SDS, 20%	_____	_____	_____

**Quality Control**

QC027 Southern Blotting and Hybridization

made by: \_\_\_\_\_ date: \_\_\_\_\_



Initials: RS

Date: 10/7/96

**S079 Hydrogen Peroxide, 3%**

lot number: \_\_\_\_\_

standard batch size: 30 X 0.5 mL

Ingredients	final concentration	amount
RM176 hydrogen peroxide, 30%	3 %	1.5 mL $\pm$ 0.1 mL
deionized water	-----	13.5 mL (guideline)

**Procedure**

Add hydrogen peroxide to a 15 mL disposable tube.

Add deionized water to a final volume of 15 mL.

Aliquot approximately 0.5 mL of hydrogen peroxide into 1.5 mL microcentrifuge tubes.

Label each tube with "H<sub>2</sub>O<sub>2</sub>" and the lot number. Label the rack with expiration date.

Store at 4°C in the dark.

Discard after 2 months.

Data Log	source	lot	amount
RM176 hydrogen peroxide, 30%	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *RSJ*

Date: *10/7/96*

**S032 Lambda Marker**

lot number: \_\_\_\_\_

page 1 of 2

Ingredients	initial concentration (ng/μL)	initial volume (μL)	final concentration	final volume (μL)
RM155 lambda Hind III fragments			20 ng/μL	
S021 yield gel loading buffer	5 X		1 X	----
S059 sterile water	-----		-----	-----

**Calculations**

Record the initial concentration in ng/μL and the initial volume in μL of the lambda Hind III DNA received from the manufacturer.

Calculate the final volume according to equation 1.

$$(\text{final volume}) = \frac{(\text{initial DNA concentration})(\text{initial DNA volume})}{(20 \text{ ng/}\mu\text{L})} \quad \text{equation 1}$$

Record the final volume above. The final volume is the total batch size.

Calculate the amount of buffer to be added according to equation 2.

$$(\text{buffer volume}) = 0.2(\text{final volume}) \quad \text{equation 2}$$

Calculate the amount of sterile water to be added according to equation 3.

$$(\text{water volume}) = [0.8 * (\text{final volume})] - (\text{initial DNA volume}) \quad \text{equation 3}$$

Record the buffer and water volumes above.

To check the calculations, add together the initial volumes of DNA, loading buffer, and sterile water.

The sum of the initial volumes must be equal to the calculated final volume.

Initials: RCJ

Date: 10/7/82

**S032 Lambda Marker**

lot number: \_\_\_\_\_

page 2 of 2

**Procedure**

Combine the DNA, loading buffer, and sterile water.

Mix well.

Using sterile pipet tips, dispense 500  $\mu$ L aliquots into sterile 1.5 mL eppendorf tubes.

Store at  $-20^{\circ}\text{C}$ .

**Data Log**

	source	lot	amount
RM155 lambda Hind III fragments	_____	_____	_____
S021 yield gel loading buffer	_____	_____	_____
S059 sterile water	_____	_____	_____

**Quality Control**

QC026 Gel Electrophoresis

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCJ

Date: 10/17/96

**S015 Lithium Chloride, 7.5 M**

lot number: \_\_\_\_\_

standard batch size: 100 mL

**Ingredients**

final  
concentration

amount

RM032 lithium chloride

7.5 M

31.8 ± 0.2 g

**Procedure**

Dissolve the lithium chloride in approximately 75 mL deionized water.

Mix well.

When the lithium chloride has dissolved, bring up to the final volume with deionized water.

Dispense 10 mL aliquots into 15 mL centrifuge tubes.

Autoclave at 250°F for 20 minutes.

Store at -20°C.

**Data Log**

source

lot

amount

RM032 lithium chloride

\_\_\_\_\_

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: PCS

Date: 10/7/86

**S055 Magnesium Chloride, 0.1M**

lot number: \_\_\_\_\_

standard batch size: 250 mL

**Ingredients**

final  
concentration

amount

RM046 magnesium chloride, hexahydrate

0.10 M

5.1 ± 0.3 g

**Procedure**

Dissolve the magnesium chloride in approximately 200 mL deionized water.

Mix well.

When the magnesium chloride has dissolved, bring up to the final volume with deionized water.

Dispense into 125 mL bottles.

Autoclave at 250°F for 20 minutes.

Store at room temperature.

**Data Log**

source

lot

amount

RM046 magnesium chloride,  
hexahydrate

\_\_\_\_\_

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: PCSDate: 10/2/96**S112 dNTP's (10mM stock)**

lot number: \_\_\_\_\_

standard batch size: ~ 25 tubes x 250  $\mu$ L**Ingredients**

	final concentration	amount
dATP, 10 mM, 320 $\mu$ L/tube	2.5 mM	1600 $\mu$ L (5 tubes)
dCTP, 10 mM, 320 $\mu$ L/tube	2.5 mM	1600 $\mu$ L (5 tubes)
dGTP, 10 mM, 320 $\mu$ L/tube	2.5 mM	1600 $\mu$ L (5 tubes)
dTTP, 10 mM, 320 $\mu$ L/tube	2.5 mM	1600 $\mu$ L (5 tubes)
Autoclaved, microcentrifuge tubes	_____	~25 tubes

**Procedure**

**NOTE: ALIQUOT ALL TUBES AT ONE TIME AND IN A ROOM FREE FROM AMPLIFIED DNA TO MINIMIZE CONTAMINATION. USING CLEAN GLOVES IS ESSENTIAL; CHANGE THEM AS OFTEN AS NEEDED. USE ONLY FILTER PIPET TIPS OR A REPEAT PIPETTOR FOR ALL PIPETTING.**

Clean the bench top thoroughly using a 10% bleach solution, and cover it with new bench paper.

Add the dNTP's together in a 15 mL sterile centrifuge tube and mix.

While wearing clean gloves, remove all tubes from the bag and place them in a clean rack designated for the PCR preparation room only.

Aliquot 250  $\mu$ L of dNTP mix into each tube.

Once aliquotting is complete, cap all tubes and store in a labeled rack away from all sources of DNA. Store at 2-8°C.

**Data Log**

	source	lot	amount
RM0211 Geneamp dNTP's	_____	_____	_____
dATP	_____	_____	_____
dCTP	_____	_____	_____
dGTP	_____	_____	_____
dTTP	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *RCJ*

Date: *10/7/96*

**S006 Phenol**

lot number: \_\_\_\_\_

standard batch size: 500 mL

page 1 of 2

Ingredients	final concentration	amount
RM112 phenol	---	500 $\pm$ 10 g
S052 TRIS, 0.1 M - pH 7.8	---	400 $\pm$ 10 mL
RM036 m-cresol	---	25 $\pm$ 1 mL
RM049 2-mercaptoethanol	---	1.0 $\pm$ 0.1 mL
RM029 4-hydroxyquinoline	---	0.50 $\pm$ 0.01 g

**Procedure**

Place a 500 g bottle of phenol crystals in a 65°C waterbath to dissolve (about 10-15 minutes).

When dissolved, add 100 mL TRIS solution to the bottle. Invert several times to mix thoroughly.

Return the bottle to the 65°C water bath and allow the temperature to equilibrate another 10-15 minutes.

Pour the solution into a separatory funnel.

When the phases have separated, return the organic phase (the bottom phase) to the bottle. Discard the aqueous phase into an organic waste container.

To the organic phase, add another 100 mL TRIS.

Separate the phases as before, isolating the organic phase (the bottom phase) each time and discarding the aqueous phase into organic waste.

Repeat the 100 mL TRIS wash one more time.

To the final organic phase, add the m-cresol, mercaptoethanol, hydroxyquinoline, and 100 mL TRIS.

Store at 2-8°C.

Initials: ACS

Date: 10/7/96

**S006 Phenol**

lot number: \_\_\_\_\_

page 2 of 2

**Data Log**

	source	lot	amount
RM112 phenol	_____	_____	_____
S052 TRIS, 0.1M - pH 7.8	_____	_____	_____
RM036 m-cresol	_____	_____	_____
RM049 2-mercaptoethanol	_____	_____	_____
RM029 4-hydroxyquinoline	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_



Initials: 29

Date: 10/7/96

**S053 Phenol-Chloroform-Isoamyl Alcohol**

lot number: \_\_\_\_\_

standard batch size: 1 L

Ingredients	final concentration	amount
S006 phenol	50 %	500 ± 5 mL
S027 chloroform-isoamyl alcohol	50 %	500 ± 5 mL

**Procedure**

**NOTE:** Use only glass containers and graduated cylinders.

Take 500 mL of the phenol mixture from the bottom phase and place it in a brown 1 L bottle.

Add 500 mL chloroform-isoamyl alcohol.

Store at 2-8°C in the flammable materials refrigerator.

Data Log	source	lot	amount
S006 phenol	_____	_____	_____
S027 chloroform-isoamyl alcohol	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *RC*

Date: *10/7/96*

**S042 Phi-X Marker**

lot number: \_\_\_\_\_

page 1 of 2

Ingredients	initial concentration (ng/μL)	initial volume (μL)	final concentration	final volume (μL)
RM156 phi-X-174, Hae III fragments			50 ng/μL	
S018 analytical gel loading buffer	5 X		1 X	----
S059 sterile water	-----		----	----

### Calculations

Record the initial concentration in ng/μL and the initial volume in μL of the phi-X-174 Hae III received from the manufacturer.

Calculate the final volume according to equation 1.

$$(\text{final volume}) = \frac{(\text{initial DNA concentration})(\text{initial DNA volume})}{(50 \text{ ng/}\mu\text{L})} \quad \text{equation 1}$$

Record the final volume above. The final volume is the total batch size.

Calculate the amount of buffer to be added according to equation 2.

$$(\text{buffer volume}) = 0.2(\text{final volume}) \quad \text{equation 2}$$

Calculate the amount of sterile water to be added according to equation 3.

$$(\text{water volume}) = [0.8 * (\text{final volume})] - (\text{initial DNA volume}) \quad \text{equation 3}$$

Record the buffer and water volumes above.

To check the calculations, add together the initial volumes of DNA, loading buffer, and sterile water.

The sum of the initial volumes must be equal to the calculated final volume.

Initials: pes

Date: 10/3/96

**S042 Phi-X Marker**

lot number: \_\_\_\_\_

page 2 of 2

### Procedure

Combine the DNA, loading buffer, and sterile water.

Mix well.

Using sterile pipet tips, dispense 500  $\mu$ L aliquots into sterile 1.5 mL eppendorf tubes.

Store at -20°C.

### Data Log

RM156 phi-X-174 Hae III fragments

S018 analytical gel loading buffer

S059 sterile water

source	lot	amount
--------	-----	--------

_____	_____	_____
-------	-------	-------

_____	_____	_____
-------	-------	-------

_____	_____	_____
-------	-------	-------

### Quality Control

QC026 Gel Electrophoresis

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCIDate: 10/3/96**S034 Phosphate Buffered Saline (PBS)**

lot number: \_\_\_\_\_

standard batch size: 4 L

Ingredients	final concentration	amount
RM005 sodium chloride	137 mM	32.0 ± 0.1 g
RM053 potassium chloride	3.0 mM	0.90 ± 0.01 g
RM065 sodium phosphate, dibasic	6.0 mM	3.41 ± 0.03 g
RM056 potassium phosphate, monobasic	1.5 mM	0.82 ± 0.02 g

**Procedure**

Add all the components to approximately 3 L deionized water.

Mix well.

Adjust the pH to 7.5.

Bring up to the final volume with deionized water.

Measure and record the final pH.

Dispense into 50 mL centrifuge tubes.

Autoclave at 250°F for 20 minutes.

Store at room temperature.

**Data Log**

	source	lot	amount
RM005 sodium chloride	_____	_____	_____
RM053 potassium chloride	_____	_____	_____
RM065 sodium phosphate, dibasic	_____	_____	_____
RM056 potassium phosphate, monobasic	_____	_____	_____

**Quality Control**

final pH: \_\_\_\_\_ spec: 7.5 ± 0.1

QC023 QuantiBlot Quality Control of Solutions- Test 150 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RC

Date: 10/17/96

**S109 PM PCR Reaction Mixture**

lot number: \_\_\_\_\_

standard batch size: ~ 65 tubes x 40  $\mu$ L

Ingredients	final concentration	amount
PM PCR reaction mix	_____	2.4 mL
Autoclaved, PCR reaction tubes	_____	55 tubes

**Procedure**

**NOTE: ALIQUOT ALL TUBES AT ONE TIME AND IN A ROOM FREE FROM AMPLIFIED DNA TO MINIMIZE CONTAMINATION. USING CLEAN GLOVES IS ESSENTIAL; CHANGE THEM AS OFTEN AS NEEDED.**

Clean the bench top thoroughly using a 10% bleach solution, and cover it with new bench paper.

While wearing clean gloves, remove all tubes from the bag and place them in a clean rack designated for the PCR preparation room only.

Using a dedicated positive displacement repeat pipettor or tips with hydrophobic filters, carefully aliquot 40  $\mu$ L of PCR reaction mixture into each tube.

Once aliquotting is complete, cap all tubes and store in a labeled rack away from all sources of DNA.

Store at 2-8°C.

Data Log	source	lot	amount
PM reaction mix	_____	_____	_____
PCR reaction tubes	_____	_____	_____

**Quality Control**

QC015 PCR Kit Amplification- Only for the first kit of each shipment/lot

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCI

Date: 10/7/96

**S126 Positive Control-External**

lot number \_\_\_\_\_

standard batch size ~50 tubes

**Procedure**

Collect EDTA blood from a volunteer. Prepare bloodstains from this blood sample.

Extract two 3x3 cm portions of the dried bloodstain following the organic extraction procedure in the RFLP manual.

Run an aliquot of each fraction on a yield gel. Pool both fractions. The DNA has to be of good quality that means of high molecular weight

Re-quantitate using the Quantiblot procedure.

Dilute with sterile water an aliquot of the DNA based on the yield gel result to a concentration of 1.25-2.5ng/20-25ul.

Quantitate dilution using Quantiblot.

<b>Data Log</b>	Source	Yield Gel Fract. 1	Yield Gel Fract. 2	Quantiblot Pool	Quantiblot Dilution
EDTA Blood	_____	_____	_____	_____	_____
	Source	Lot	Amount		
S059 Sterile Water	_____	_____	_____		

**Quality Control**

QC031- QUAD STR/PCR Amplification

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RC

Date: 10/7/96

## **S127 Positive Control-Internal**

batch number \_\_\_\_\_

**page 1 of 2**

standard batch size 24 tubes

### **Procedure**

All previously typed samples can be used as an internal positive control (reanalysis sample) if they meet the following criteria:

- the sample has been successfully typed for all STR loci so that all alleles for each locus are known.
- the DNA concentration according to Quantiblot was greater or equal to 1.25ng/20ul.
- a sufficient amount of the DNA extract is left.
- the peak height for VWA was greater or equal 200 f.u. during the first analysis.

The internal positive control samples are prepared by **two** analysts together as follows:

Select 24 appropriate DNA samples

Fill out the Positive Control-Internal solution sheet. Note: One person should read the allele designations from the electropherograms, not the table, while the second person fills out the sheet

Calculate the amounts of DNA and deionized water according to the STR Quad amplification sheet.

Label 24 tubes adhering to this format: 1.1 to 1.24 for the first batch, 2.1 to 2.24 for the second batch and so on. Label tubes for 1:10 dilutions if required.

Pipet the calculated amount of water in all tubes and close the tubes.

Spin down Chelex extracts. With one analyst watching, pipet the required amount of DNA solution into the correct tube. Make 1:10 dilutions and use these where required.

Mix tubes, spin down shortly, and place them in the Positive Control Internal rack in the PCR set up room.

Initials: RS

Date: 10/7/96

**S127 Positive Control-Internal**

batch number \_\_\_\_\_

**page 2 of 2**

Batch. PI #	MB / ex. #	Conc. ng/20 $\mu$ l	DNA vol $\mu$ l	dH <sub>2</sub> O vol $\mu$ l	VWA	F13A1	THO1	FES
.1								
.2								
.3								
.4								
.5								
.6								
.7								
.8								
.9								
.10								
.11								
.12								
.13								
.14								
.15								
.16								
.17								
.18								
.19								
.20								
.21								
.22								
.23								
.24								

S059 Sterile Water \_\_\_\_\_

Prepared by \_\_\_\_\_ Witness \_\_\_\_\_ Date \_\_\_\_\_

October 3, 1996



Initials: RCJ

Date: 10/7/96

**S097 Pre-Wetting Solution**

lot number: \_\_\_\_\_

standard batch size: 4 L

Ingredients	final concentration	amount
RM004 NaOH, 10 N	0.4 N	160 ± 10 mL
S009 EDTA, 0.5 M	25 mM	200 ± 10 mL

**Procedure**

Measure 3640 mL deionized water into a 4 L bottle.

Add 160 mL NaOH and 200 mL EDTA.

Cap and mix well by inverting.

Dispense into 1 L bottles or store in bulk.

Store at room temperature.

Data Log	source	lot	amount
RM004 NaOH, 10 N	_____	_____	_____
S009 EDTA, 0.5 M	_____	_____	_____

made by: \_\_\_\_\_

date: \_\_\_\_\_

Initials: 24

Date: 10/7/96

**S118 Primer, VWA 1 (50 µM)**

lot number: \_\_\_\_\_

**page 1 of 2**

**Physical data**

Sequence JOE - 5' CC CTA GTG GAT GAT AAG AAT AAT CAG TAT 3'

Oligo	M.W.	µg/ O.D.	pmol/ O.D.
VWA 1	9272.0	30.1	3246.3

Ingredients	initial amount (O.D.)	amount in pmoles	final concentration	volume H <sub>2</sub> O (µL)
RM 1087 VWA 1 primer			50 pM/µL	-----
S059 Sterile Water	-----	-----	-----	

**Calculations**

Record the initial amount in O.D. received from the manufacturer.

Calculate the total amount in pmoles according to equation 1.

$$(\text{Amount in pmoles}) = (\text{O.D.}) \times 3246 \quad \text{equation 1}$$

Record the amount in pmoles above.

Calculate the amount of dH<sub>2</sub>O to be added according to equation 2.

$$(\text{dH}_2\text{O volume}) = \frac{(\text{amount in pmoles})}{50} \quad \text{equation 2}$$

Record the water volume above.

Have somebody check the calculation.

Initials: RC

Date: 10/7/96

**S118 Primer, VWA 1 (50 µM)**

lot number: \_\_\_\_\_  
**page 2 of 2**

**Procedure**

Add the sterile water to the original primer tube.

Mix well.

Dispense 50 µL aliquots into 1.8 mL microcentrifuge tubes.

Store at -20°C.

**Data Log**

	source	lot	amount
RM1087 Primer VWA 1	_____	_____	_____
S059 sterile water	_____	_____	_____

**Calculation checked by** \_\_\_\_\_

**Quality Control**

QC031 QUAD STR/PCR Amplification

QC023 Quantiblot- test 1µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *Ad*

Date: *10/7/96*

**S119 Primer, VWA 2 (50 µM)**

lot number: \_\_\_\_\_  
**page 1 of 2**

**Physical data**

Sequence 5' GGA CAG ATG ATA AAT ACA TAG GAT GGA TGG 3'

Oligo	M.W.	µg/ O.D.	pmol/ O.D.
VWA 2	9383.0	29.4	3133.3

Ingredients	initial amount (O.D.)	amount in pmoles	final concentration	volume H <sub>2</sub> O (µL)
RM 1080 VWA 2 primer			50 pM/µL	-----
S059 Sterile Water	-----	-----	-----	

**Calculations**

Record the initial amount in O.D. received from the manufacturer.

Calculate the total amount in pmoles according to equation 1.

$$(\text{Amount in pmoles}) = (\text{O.D.}) \times 3133 \quad \text{equation 1}$$

Record the amount in pmoles above.

Calculate the amount of dH<sub>2</sub>O to be added according to equation 2.

$$(\text{dH}_2\text{O volume}) = \frac{(\text{amount in pmoles})}{50} \quad \text{equation 2}$$

Record the water volume above.

Have somebody check the calculation.

Initials: RC

Date: 10/7/96

**S119 Primer, VWA 2 (50  $\mu$ M)**

lot number: \_\_\_\_\_  
**page 2 of 2**

**Procedure**

Add the sterile water to the original primer tube.

Mix well.

Dispense 50  $\mu$ L aliquots into 1.8 mL microcentrifuge tubes.

Store at -20°C.

**Data Log**

source	lot	amount
--------	-----	--------

RM1080 Primer VWA 2	_____	_____
---------------------	-------	-------

S059 sterile water	_____	_____
--------------------	-------	-------

**Calculation checked by** \_\_\_\_\_

**Quality Control**

QC031 QUAD STR/PCR Amplification

QC023 Quantiblot- test 1 $\mu$ L of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *ACJ*

Date: *10/7/96*

**S120 Primer THO1 1 (50 µM)**

lot number: \_\_\_\_\_

**page 1 of 2**

**Physical data**

Sequence FAM - 5' GT GGG CTG AAA AGC TCC CGA TTA T 3'

Oligo	M.W.	µg/ O.D.	pmol/ O.D.
THO1 1	7386.1	32.3	4373.1

Ingredients	initial amount (O.D.)	amount in pmoles	final concentration	volume H <sub>2</sub> O (µL)
RM 1086 Tho1 1 primer			50 pM/µL	-----
S059 Sterile Water	-----	-----	-----	

**Calculations**

Record the initial amount in O.D. received from the manufacturer.

Calculate the total amount in pmoles according to equation 1.

$$(\text{Amount in pmoles}) = (\text{O.D.}) \times 4373 \quad \text{equation 1}$$

Record the amount in pmoles above.

Calculate the amount of dH<sub>2</sub>O to be added according to equation 2.

$$(\text{dH}_2\text{O volume}) = \frac{(\text{amount in pmoles})}{50} \quad \text{equation 2}$$

Record the water volume above.

Have somebody check the calculation.

Initials: CS

Date: 10/7/96

**S120 Primer, Tho1 1 (50  $\mu$ M)**

lot number: \_\_\_\_\_  
**page 2 of 2**

### Procedure

Add the sterile water to the original primer tube.

Mix well.

Dispense 50  $\mu$ L aliquots into 1.8 mL microcentrifuge tubes.

Store at -20°C.

### Data Log

source	lot	amount
--------	-----	--------

RM1086 Primer Tho1 1	_____	_____
----------------------	-------	-------

S059 sterile water	_____	_____
--------------------	-------	-------

**Calculation checked by** \_\_\_\_\_

### Quality Control

QC031 QUAD STR/PCR Amplification

QC023 Quantiblot- test 1  $\mu$ L of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *RV*

Date: *10/7/96*

**S121 Primer THO1 2 (50  $\mu$ M)**

lot number: \_\_\_\_\_  
page 1 of 2

**Physical data**

Sequence 5' GTG ATT CCC ATT GGC CTG TTC CTC 3'

Oligo	M.W.	$\mu$ g/ O.D.	pmol/ O.D.
THO1 2	7257.8	35.1	4836.2

Ingredients	initial amount (O.D.)	amount in pmoles	final concentration	volume H <sub>2</sub> O ( $\mu$ L)
RM 1079Tho1 2 primer			50 pM/ $\mu$ L	-----
S059 Sterile Water	-----	-----	-----	

**Calculations**

Record the initial amount in O.D. received from the manufacturer.

Calculate the total amount in pmoles according to equation 1.

$$(\text{Amount in pmoles}) = (\text{O.D.}) \times 4836 \quad \text{equation 1}$$

Record the amount in pmoles above.

Calculate the amount of dH<sub>2</sub>O to be added according to equation 2.

$$(\text{dH}_2\text{O volume}) = \frac{(\text{amount in pmoles})}{50} \quad \text{equation 2}$$

Record the water volume above.

Have somebody check the calculation.



Initials: RC

Date: 10/7/96

**S121 Primer, Tho1 2 (50 µM)**

lot number: \_\_\_\_\_  
**page 2 of 2**

**Procedure**

Add the sterile water to the original primer tube.

Mix well.

Dispense 50 µL aliquots into 1.8 mL microcentrifuge tubes.

Store at -20°C.

**Data Log**

source	lot	amount
--------	-----	--------

RM1079 Primer Tho1 2	_____	_____
----------------------	-------	-------

S059 sterile water	_____	_____
--------------------	-------	-------

**Calculation checked by** \_\_\_\_\_

**Quality Control**

QC031 QUAD STR/PCR Amplification

QC023 Quantiblot- test 1µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *RC*

Date: *10/7/96*

**S122 Primer F13A1 1 (50  $\mu$ M)**

lot number: \_\_\_\_\_

**page 1 of 2**

**Physical data**

Sequence JOE - 5' AT GCC ATG CAG ATT AGA AA 3'

Oligo	M.W.	$\mu$ g/ O.D.	pmol/ O.D.
F13A1/1	5841.8	29.8	5101.2

Ingredients	initial amount (O.D.)	amount in pmoles	final concentration	volume H <sub>2</sub> O ( $\mu$ L)
RM 1085 F13A1 1 primer			50 pM/ $\mu$ L	-----
S059 Sterile Water	-----	-----	-----	

**Calculations**

Record the initial amount in O.D. received from the manufacturer.

Calculate the total amount in pmoles according to equation 1.

$$(\text{Amount in pmoles}) = (\text{O.D.}) \times 5101 \quad \text{equation 1}$$

Record the amount in pmoles above.

Calculate the amount of dH<sub>2</sub>O to be added according to equation 2.

$$(\text{dH}_2\text{O volume}) = \frac{(\text{amount in pmoles})}{50} \quad \text{equation 2}$$

Record the water volume above.

Have somebody check the calculation.

Initials: AS

Date: 10/7/96

**S122 Primer, F13A1 1 (50  $\mu$ M)**

lot number: \_\_\_\_\_

**page 2 of 2**

### Procedure

Add the sterile water to the original primer tube.

Mix well.

Dispense 50  $\mu$ L aliquots into 1.8 mL microcentrifuge tubes.

Store at -20°C.

### Data Log

source

lot

amount

RM1085 Primer, F13A1 1

\_\_\_\_\_

S059 sterile water

\_\_\_\_\_

**Calculation checked by** \_\_\_\_\_

### Quality Control

QC031 QUAD STR/PCR Amplification

QC023 Quantiblot- test 1  $\mu$ L of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *24*

Date: *06/2/96*

**S123 Primer, F13A1 2 (50 µM)**

lot number: \_\_\_\_\_

**page 1 of 2**

**Physical data**

Sequence 5' GAG GTT GCA CTC CAG CCT TT 3'

Oligo	M.W.	µg/ O.D.	pmol/ O.D.
F13A1/2	6080.0	34.1	5608.6

Ingredients	initial amount (O.D.)	amount in pmoles	final concentration	volume H <sub>2</sub> O (µL)
RM 1076 F13A1 2 primer			50 pM/µL	-----
S059 Sterile Water	-----	-----	-----	

**Calculations**

Record the initial amount in O.D. received from the manufacturer.

Calculate the total amount in pmoles according to equation 1.

$$(\text{Amount in pmoles}) = (\text{O.D.}) \times 5609 \quad \text{equation 1}$$

Record the amount in pmoles above.

Calculate the amount of dH<sub>2</sub>O to be added according to equation 2.

$$(\text{dH}_2\text{O volume}) = \frac{(\text{amount in pmoles})}{50} \quad \text{equation 2}$$

Record the water volume above.

Have somebody check the calculation.

Initials: RCJ

Date: 10/7/96

**S123 Primer, F13A1 2 (50  $\mu$ M)**

lot number: \_\_\_\_\_

**page 2 of 2**

### **Procedure**

Add the sterile water to the original primer tube.

Mix well.

Dispense 50  $\mu$ L aliquots into 1.8 mL microcentrifuge tubes.

Store at -20°C.

### **Data Log**

source

lot

amount

RM 1076 Primer F13A1 2

\_\_\_\_\_

S059 sterile water

\_\_\_\_\_

**Calculation checked by** \_\_\_\_\_

### **Quality Control**

QC031 QUAD STR/PCR Amplification

QC023 Quantiblot- test 1 $\mu$ L of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *RCJ*

Date: *10/7/96*

**S124 Primer, FES/FPS 1 (50 µM)**

lot number: \_\_\_\_\_  
**page 1 of 2**

**Physical data**

Sequence 5' GG GAT TTC CCT ATG GAT TGG 3'

Oligo	M.W.	µg/ O.D.	pmol/ O.D.
FES 1	6173	32.8	5313.5

Ingredients	initial amount (O.D.)	amount in pmoles	final concentration	volume H <sub>2</sub> O (µL)
RM 1084 FES 1 primer			50 pM/µL	-----
S059 Sterile Water	-----	-----	-----	

**Calculations**

Record the initial amount in O.D. received from the manufacturer.

Calculate the total amount in pmoles according to equation 1.

$$(\text{Amount in pmoles}) = (\text{O.D.}) \times 5314 \quad \text{equation 1}$$

Record the amount in pmoles above.

Calculate the amount of dH<sub>2</sub>O to be added according to equation 2.

$$(\text{dH}_2\text{O volume}) = \frac{(\text{amount in pmoles})}{50} \quad \text{equation 2}$$

Record the water volume above.

Have somebody check the calculation.

Initials: RCS

Date: 10/7/96

**S124 Primer, FES/FPS (50  $\mu$ M)**

lot number: \_\_\_\_\_

**page 2 of 2**

### **Procedure**

Add the sterile water to the original primer tube.

Mix well.

Dispense 50  $\mu$ L aliquots into 1.8 mL microcentrifuge tubes.

Store at -20°C.

### **Data Log**

source

lot

amount

RM1084 Primer FES/FPS 1

\_\_\_\_\_

S059 sterile water

\_\_\_\_\_

**Calculation checked by** \_\_\_\_\_

### **Quality Control**

QC031 QUAD STR/PCR Amplification

QC023 Quantiblot- test 1  $\mu$ L of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *FCS*

Date: *10/7/86*

**S125 Primer FES/FPS 2 (50  $\mu$ M)**

lot number: \_\_\_\_\_  
**page 1 of 2**

**Physical data**

Sequence FAM - 5' GCG AAA GAA TGA GAC TAC AT 3'

Oligo	M.W.	$\mu$ g/ O.D.	pmol/ O.D.
FES 2	6179	29.7	4806.6

Ingredients	initial amount (O.D.)	amount in pmoles	final concentration	volume H <sub>2</sub> O ( $\mu$ L)
RM 1075 FES 2 primer			50 pM/ $\mu$ L	-----
S059 Sterile Water	-----	-----	-----	

**Calculations**

Record the initial amount in O.D. received from the manufacturer.

Calculate the total amount in pmoles according to equation 1.

$$(\text{Amount in pmoles}) = (\text{O.D.}) \times 4807 \quad \text{equation 1}$$

Record the amount in pmoles above.

Calculate the amount of dH<sub>2</sub>O to be added according to equation 2.

$$(\text{dH}_2\text{O volume}) = \frac{(\text{amount in pmoles})}{50} \quad \text{equation 2}$$

Record the water volume above.

Have somebody check the calculation.



Initials: RCJ

Date: 10/7/96

**S125 Primer, FES/FPS (50  $\mu$ M)**

lot number: \_\_\_\_\_

page 2 of 2

### Procedure

Add the sterile water to the original primer tube.

Mix well.

Dispense 50  $\mu$ L aliquots into 1.8 mL microcentrifuge tubes.

Store at -20°C.

### Data Log

source

lot

amount

RM1075 Primer FES/FPS 2

\_\_\_\_\_

S059 sterile water

\_\_\_\_\_

Calculation checked by \_\_\_\_\_

### Quality Control

QC031 QUAD STR/PCR Amplification

QC023 Quantiblot- test 1  $\mu$ L of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCJ

Date: 10/7/96

**S011 Protein Lysis Buffer (PLB)**

lot number: \_\_\_\_\_

standard batch size: 2 L

Ingredients	final concentration	amount
S009 EDTA, 0.5M	10 mM	40 ± 2 mL
RM073 TRIS	10 mM	2.4 ± 0.1 g
RM005 sodium chloride	10 mM	1.2 ± 0.05 g
RM096 hydrochloric acid	---	

**Procedure**

Add the TRIS, EDTA, and sodium chloride to approximately 1.5 L deionized water.  
Mix well.

Adjust the pH to 7.9 with hydrochloric acid

Raise to the final volume with deionized water.

Mix well.

Dispense into 15 mL centrifuge tubes.

Autoclave at 250°F for 30 minutes.

Store at 2-8°C.

Data Log	source	lot	amount
S009 EDTA, 0.5M	_____	_____	_____
RM073 TRIS	_____	_____	_____
RM005 sodium chloride	_____	_____	_____
RM096 hydrochloric acid	_____	_____	_____

**Quality Control**

final pH: \_\_\_\_\_ spec: 7.9 ± 0.1

QC023 QuantiBlot Quality Control of Solutions- test 150 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RSJ

Date: 10/7/96

**S014 Proteinase-K Enzyme, 10mg/mL**

lot number: \_\_\_\_\_

standard batch size: 10 mL

**Ingredients**

final  
concentration

amount

RM119 proteinase-K, lyophilized

10 mg/mL

100 ± 1 mg

S059 sterile water

\_\_\_\_\_

10 ± 0.5 mL

**Procedure**

Add 10 mL sterile, deionized water to one bottle (100 mg) lyophilized proteinase-K enzyme.

Mix by slowly inverting until completely reconstituted.

Dispense 500 µL aliquots into 1.5 mL eppendorf tubes.

Store at -20°C.

**Data Log**

source

lot

amount

RM119 proteinase-K, lyophilized

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

S059 sterile water

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**Quality Control**

QC023 QuantiBlot Quality Control of Solutions- Test 10 µL of solution

QC024 Non-Organic Extraction

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCJ

Date: 10/7/96

**S037 Proteinase-K Enzyme, 20mg/mL**

lot number: \_\_\_\_\_

standard batch size: 5 mL

**Ingredients**

final  
concentration

amount

RM119 proteinase-k, lyophilized

20 mg/mL

100.0 ± 0.5 mg

S059 sterile water

—

5 mL

**Procedure**

Add 5 mL sterile water to one bottle (100 mg) lyophilized proteinase-k enzyme.

Mix by slowly inverting until completely reconstituted.

Dispense 500 ul aliquots into 1.5 mL eppendorf tubes.

Store at -20°C.

**Data Log**

source

lot

amount

RM119 proteinase-k, lyophilized

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

S059 sterile water

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**Quality Control**

QC023 QuantiBlot Quality Control of Solutions- Test 10 µL of solution

QC024 Non-Organic Extraction

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *QC*

Date: *10/7/96*

**S111 QUAD STR/PCR Reaction Mixture**

standard batch size: 50-200 tubes

lot number: \_\_\_\_\_

**page 1 of 2**

<b><u>Ingredients:</u></b>	<b><u>Final</u></b>	<b><u>1 Tube</u></b>	<b><u>50</u></b>	<b><u>100</u></b>	<b><u>200</u></b>
	<b><u>Concentration</u></b>	<b><u>Amount</u></b>	<b><u>Tubes</u></b>	<b><u>Tubes</u></b>	<b><u>Tubes</u></b>
RM0275 10X PCR Buffer II	1X	5 µL	250 µL	500 µL	1000 µL
S112 dNTP's (2.5 mM)	200 mM	4 µL	200 µL	400 µL	800 µL
S059 sterile dH2O	-----	6.8 µL	340 µL	680 µL	1360 µL
S117 BSA (5mg/mL)	160ug/ml	1.6 µL	80 µL	160 µL	320 µL
S118 VWA/1 (50pM/µL)	0.22 mM	0.22 µL	11 µL	22 µL	44 µL
S119 VWA/2 (50pM/µL)	0.22 mM	0.22 µL	11 µL	22 µL	44 µL
S120 THO1/1 (50pM/µL)	0.22 mM	0.22 µL	11 µL	22 µL	44 µL
S121 THO1/2 (50pM/µL)	0.22 mM	0.22 µL	11 µL	22 µL	44 µL
S122 F13A1/1 (50pM/µL)	0.16 mM	0.16 µL	8 µL	16 µL	32 µL
S123 F13A1/2 (50pM/µL)	0.16 mM	0.16 µL	8 µL	16 µL	32 µL
S124 FES/1/(50pM/µL)	0.20 mM	0.20 µL	10 µL	20 µL	40 µL
S125 FES/2 (50pM/µL)	0.20 mM	0.20 µL	10 µL	20 µL	40 µL
<u>RM0275 AmpliTaq (5u/µL)</u>	5 U	<u>1 µL</u>	<u>50 µL</u>	<u>100 µL</u>	<u>200 µL</u>
<b>TOTAL</b>		<b>20 µL</b>	<b>1 mL</b>	<b>2 mL</b>	<b>4 mL</b>

**Procedure**

**NOTE: ALIQUOT ALL TUBES AT ONE TIME AND IN A ROOM FREE FROM AMPLIFIED DNA TO MINIMIZE CONTAMINATION. USING CLEAN GLOVES IS ESSENTIAL; CHANGE THEM AS OFTEN AS NEEDED.**

Clean the bench top thoroughly using a 10% bleach solution, and cover it with new bench paper. Add the ingredients to either a microcentrifuge tube or a 15 mL centrifuge tube using pipetmen dedicated to PCR preparation area only.

While wearing clean gloves, remove sufficient amount of tubes from the bag and place them in a clean rack designated for the PCR prep room only.

Vortex and spin briefly. Add 20 µL per tube using a dedicated repeat pipettor or tips with hydrophobic filters.

Cap all tubes and store in a labeled rack away from all sources of DNA.

Store at 2-8°C.

Initials: *RG*

Date: *10/7/96*

**S111 QUAD STR/PCR Reaction Mixture**

lot number: \_\_\_\_\_

**page 2 of 2**

**Data Log**

	source	lot	amount
RM0275 10X PCR Buffer II	_____	_____	_____
S112 dNTP's (2.5 mM)	_____	_____	_____
S059 sterile dH2O	_____	_____	_____
S117 BSA (5mg/mL)	_____	_____	_____
S118 VWA/1 (50pM/μL)	_____	_____	_____
S119 VWA/2 (50pM/μL)	_____	_____	_____
S120 THO1/1 (50pM/μL)	_____	_____	_____
S121 THO1/2 (50pM/μL)	_____	_____	_____
S122 F13A1/1 (50pM/μL)	_____	_____	_____
S123 F13A1/2 (50pM/μL)	_____	_____	_____
S124 FES/1/(50pM/μL)	_____	_____	_____
S125 FES/2 (50pM/μL)	_____	_____	_____
RM0275 AmpliTaq (5u/μL)	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *RCJ*

Date: *10/24/86*

### **S100 QuantiBlot DNA Standards**

standard batch size: variable

lot number: \_\_\_\_\_  
page 1 of 2

<b>Ingredients</b>	<b>final concentration</b>	<b>amount</b>
RM442 DNA Standard A	varies	varies
S128 TE <sup>-4</sup> , 1X	1X	varies

### **Procedure**

Each lot of QuantiBlot DNA Standards is prepared by pooling up to 10 DNA Standard A's (from the QuantiBlot kit) and serially diluting according to the following procedure:

1. Pool the contents of five or ten DNA Standard A tubes (use all one lot number).
2. Vortex to mix thoroughly.
3. Label seven sterile microfuge tubes, 1A - 1G.
4. If five DNA Standard A tubes were pooled:  
Transfer 600  $\mu$ L of DNA Standard A into the tube labeled 1A. This is now DNA Standard 1A.  
Aliquot 300  $\mu$ L of 1X TE<sup>-4</sup> into each of the six remaining tubes labeled 1B-1G.  
Add 300  $\mu$ L of DNA Standard 1A to the 300  $\mu$ L of 1X TE<sup>-4</sup> in tube 1B. Vortex to mix thoroughly.  
Add 300  $\mu$ L of diluted DNA Standard (tube 1B) to the 300  $\mu$ L of 1X TE<sup>-4</sup> in tube 1C. Vortex to mix thoroughly.  
Add 300  $\mu$ L of diluted DNA Standard (tube 1C) to the 300  $\mu$ L of 1X TE<sup>-4</sup> in tube 1D. Vortex to mix thoroughly.  
Continue the serial dilution through tube 1G.
5. If ten DNA Standard A tubes were pooled:  
Transfer 1200  $\mu$ L of DNA Standard A into the tube labeled 1A. This is now DNA Standard 1A.  
Aliquot 600  $\mu$ L of 1X TE<sup>-4</sup> into each of the six remaining tubes labeled 1B-1G.  
Add 600  $\mu$ L of DNA Standard 1A to the 600  $\mu$ L of 1X TE<sup>-4</sup> in tube 1B. Vortex to mix thoroughly.

Initials: RCJ

Date: 10/21/96

**S100 QuantiBlot DNA Standards**

standard batch size: variable

lot number: \_\_\_\_\_  
page 2 of 2

Add 600  $\mu$ L of diluted DNA Standard (tube 1B) to the 600  $\mu$ L of 1X TE<sup>-4</sup> in tube 1C.  
Vortex to mix thoroughly.

Add 600  $\mu$ L of diluted DNA Standard (tube 1C) to the 600  $\mu$ L of 1X TE<sup>-4</sup> in tube 1D.  
Vortex to mix thoroughly.

Continue the serial dilution through tube 1G.

6. Store at 2° to 8°C.
7. DNA Standards are stable for at least 3 months as 2° to 8°C.

If the dilution steps are performed as described above, the seven DNA Standard tubes will have the following concentrations of human DNA:

DNA Standards		
Standard Tube	Conc (ng/ $\mu$ L)	Quantity (ng/5 $\mu$ L)
1A	2	10
1B	1	5
1C	0.5	2.5
1D	0.25	1.25
1E	0.125	0.625
1F	0.0625	0.3125
1G	0.03125	0.15625

**Data Log**

source                      lot                      amount

RM442 DNA Standard A                      \_\_\_\_\_

S128 TE<sup>-4</sup>, 1X                      \_\_\_\_\_

**Quality Control**

QC018 QuantiBlot Hybridization.

made by: \_\_\_\_\_

date: \_\_\_\_\_



Initials: *RG*

Date: *10/17/96*

**S100 QuantiBlot DNA Standards**

standard batch size: variable

lot number: \_\_\_\_\_

page 1 of 2

Ingredients	final concentration	amount
RM442 DNA Standard A	varies	varies
S039 TE, 1X	1X	varies

**Procedure**

Each lot of QuantiBlot DNA Standards is prepared by pooling up to 10 DNA Standard A's (from the QuantiBlot kit) and serially diluting according to the following procedure:

1. Pool the contents of five or ten DNA Standard A tubes (use all one lot number).
2. Vortex to mix thoroughly.
3. Label seven sterile microfuge tubes, 1A - 1G.
4. If five DNA Standard A tubes were pooled:  
Transfer 600  $\mu$ L of DNA Standard A into the tube labeled 1A. This is now DNA Standard 1A.  
Aliquot 300  $\mu$ L of 1X TE into each of the six remaining tubes labeled 1B-1G.  
Add 300  $\mu$ L of DNA Standard 1A to the 300  $\mu$ L of 1X TE in tube 1B. Vortex to mix thoroughly.  
Add 300  $\mu$ L of diluted DNA Standard (tube 1B) to the 300  $\mu$ L of 1X TE in tube 1C. Vortex to mix thoroughly.  
Add 300  $\mu$ L of diluted DNA Standard (tube 1C) to the 300  $\mu$ L of 1X TE in tube 1D. Vortex to mix thoroughly.  
Continue the serial dilution through tube 1G.
5. If ten DNA Standard A tubes were pooled:  
Transfer 1200  $\mu$ L of DNA Standard A into the tube labeled 1A. This is now DNA Standard 1A.  
Aliquot 600  $\mu$ L of 1X TE into each of the six remaining tubes labeled 1B-1G.  
Add 600  $\mu$ L of DNA Standard 1A to the 600  $\mu$ L of 1X TE in tube 1B. Vortex to mix thoroughly.

Initials: RS

Date: 10/7/96

**S100 QuantiBlot DNA Standards**

standard batch size: variable

lot number: \_\_\_\_\_  
page 2 of 2

Add 600  $\mu$ L of diluted DNA Standard (tube 1B) to the 600  $\mu$ L of 1X TE in tube 1C.  
Vortex to mix thoroughly.

Add 600  $\mu$ L of diluted DNA Standard (tube 1C) to the 600  $\mu$ L of 1X TE in tube 1D.  
Vortex to mix thoroughly.

Continue the serial dilution through tube 1G.

6. Store at 2° to 8°C.

7. DNA Standards are stable for at least 3 months as 2° to 8°C.

If the dilution steps are performed as described above, the seven DNA Standard tubes will have the following concentrations of human DNA:

DNA Standards		
Standard Tube	Conc (ng/ $\mu$ L)	Quantity (ng/5 $\mu$ L)
1A	2	10
1B	1	5
1C	0.5	2.5
1D	0.25	1.25
1E	0.125	0.625
1F	0.0625	0.3125
1G	0.03125	0.15625

**Data Log**

source                      lot                      amount

RM442 DNA Standard A

\_\_\_\_\_

S039 TE, 1X

\_\_\_\_\_

**Quality Control**

QC018 QuantiBlot Hybridization.

made by: \_\_\_\_\_

date: \_\_\_\_\_

Initials: RS

Date: 10/7/96

**S099 QuantiBlot Wash Solution**

lot number: \_\_\_\_\_

standard batch size: 4 x 4 L

Ingredients	final concentration	amount/ 4 Liter
S002 SSPE, 20X	1.5 X	300 ± 10 mL
S001 SDS, 20%	0.5 %	100 ± 5 mL

**Procedure**

Measure 3600 mL deionized water into four 4 L bottles.

Add 300 mL SSPE and 100 mL SDS to each bottle.

Cap and mix well by inverting.

Store at room temperature.

Data Log	source	lot	amount
S002 SSPE, 20X	_____	_____	_____
S001 SDS, 20%	_____	_____	_____

made by: \_\_\_\_\_

date: \_\_\_\_\_

Initials: RCJ

Date: 10/7/96

**S054 Sarkosyl, 10%**

lot number: \_\_\_\_\_

standard batch size: 100 mL

**Ingredients**

final  
concentration

amount

S040 sarkosyl, 20%

10. %

50 ± 2 mL

**Procedure**

Dilute 50 mL of 20% sarkosyl with 50 mL deionized water.

Mix well.

Filter sterilize.

Dispense into sterile 15 mL tubes.

Store at 2-8°C.

**Data Log**

source

lot

amount

S040 sarkosyl, 20%

\_\_\_\_\_

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RC

Date: 10/7/86

**S040 Sarkosyl, 20%**

lot number: \_\_\_\_\_

standard batch size: 100 mL

Ingredients	final concentration	amount
RM057 sarkosyl	20 %	20 ± 0.5 g

**Procedure**

Add the sarkosyl to approximately 75 mL deionized water.

Mix until the solution is completely clear.

Bring up to volume with deionized water.

Filter sterilize.

Dispense into sterile 15 mL tubes.

Store at 2-8°C.

Data Log	source	lot	amount
----------	--------	-----	--------

RM057 sarkosyl	_____	_____	_____
----------------	-------	-------	-------

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: PCJ

Date: 2/9/98

**S132 Sequencing loading buffer**

lot number: \_\_\_\_\_

standard batch size: 25mL

Ingredients	final	amount
S009 500mM EDTA, pH8.0	25 mM	1.25 $\pm$ 0.05mL
RM Blue Dextran	50 mg/mL	1250 mg $\pm$ 10 mg

**Procedure**

**NOTE: PREPARE AWAY FROM AMPLIFIED DNA TO MINIMIZE CONTAMINATION. USING CLEAN GLOVES IS ESSENTIAL; CHANGE THEM AS OFTEN AS NEEDED.**

Clean the bench top thoroughly using a 10% bleach solution, and cover it with new bench paper.

Pipette EDTA into a 25mL cylinder. Fill up to 25mL using deionized water.

Decant into an 100mL Erlenmeyer flask. Add Blue Dextran. Stir at room temperature until dissolved.

Label 25 1.5mL reaction tubes.

Add 1000 $\mu$ L of the sequencing loading buffer to each tube. Close all tubes.

Store at 2-8°C.

**Data Log**

	source	lot	amount
S009 500mM EDTA, pH8.0	_____	_____	_____
RM Blue Dextran	_____	_____	_____

**Quality Control**

QC033 STR gel electrophoresis

made by: \_\_\_\_\_ date: \_\_\_\_\_

February 6, 1998

79 A

Initials: RSJ

Date: 10/7/96

**S101 SDS, 0.1%**

lot number: \_\_\_\_\_

standard batch size: 20 L

<b>Ingredients</b>	<b>final concentration</b>	<b>amount</b>
S001 SDS, 20%	0.1 %	100 ± 10 mL

### **Procedure**

Add approximately 15 L of deionized water into a 20 L carboy.

Add 100 mL 20% SDS.

Mix .

Bring up to a final volume of 20 L with deionized water.

Mix.

Store at room temperature.

<b>Data Log</b>	<b>source</b>	<b>lot</b>	<b>amount</b>
S001 SDS, 20%	_____	_____	_____

made by: \_\_\_\_\_

date: \_\_\_\_\_

Initials: RCJ

Date: 10/7/96

**S045 SDS, 10%**

standard batch size: 100 mL

lot number: \_\_\_\_\_

**Ingredients**

	final concentration	amount
RM007 sodium dodecyl sulfate	10 %	10.0 ± 0.3 g

**Procedure**

CAUTION: AN AEROSOL MASK OR FUME HOOD MUST BE USED WHEN MAKING THIS SOLUTION.. WEAR GOGGLES FOR EYE PROTECTION.

Dissolve the SDS in approximately 75 mL deionized water.

Warm the solution until all the solids have dissolved and the solution is clear.

Bring up to volume with deionized water.

Filter sterilize the warm solution.

Dispense into sterile 100 mL bottles.

Store at room temperature.

**Data Log**

	source	lot	amount
RM007 sodium dodecyl sulfate	_____	_____	_____

**Quality Control**

QC023 QuantiBlot Quality Control of Solutions- Test 25 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_



Initials: RCJ

Date: 10/7/86

**S001 SDS, 20%**

standard batch size: 1 L

lot number: \_\_\_\_\_

**Ingredients**

final  
concentration

amount

RM007 sodium dodecyl sulfate

20 %

200 ± 5 g

**Procedure**

CAUTION: AN AEROSOL MASK OR FUME HOOD MUST BE USED WHEN MAKING THIS SOLUTION. WEAR GOGGLES FOR EYE PROTECTION.

Warm approximately 750 mL deionized water on a stirring hot plate.

Add a fraction of the SDS, allowing the solids to dissolve before adding more.

Add the SDS until it is all in solution.

When the solution is clear, bring up to volume with deionized water.

Filter sterilize the warm solution.

Dispense into sterile 500 mL bottles.

Store at room temperature.

**Data Log**

source

lot

amount

RM007 SDS

\_\_\_\_\_

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCJ

Date: 10/7/85

**S046 SLS, 20%**

standard batch size: 4 L

lot number: \_\_\_\_\_

**Ingredients**

	final concentration	amount
RM218 sodium lauryl sulfate	20 %	800 ± 5 g

**Procedure**

CAUTION: AN AEROSOL MASK OR FUME HOOD MUST BE USED WHEN MAKING THIS SOLUTION. WEAR GOGGLES FOR EYE PROTECTION.

Warm approximately 2.5 L deionized water on a stirring hot plate.

Add a fraction of the SLS, allowing the solids to dissolve before adding more.

Add SLS until it is all in solution.

Mix well.

When the solution is clear, bring up to volume with deionized water.

Sterile filter each 2 L with a sterile unexpired cellulose nitrate filter.

Dispense into 1 L bottles.

Store at room temperature.

**Data Log**

	source	lot	amount
RM218 sodium lauryl sulfate	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RC

Date: 10/7/96

**S043 Sodium Acetate, 2 M**  
standard batch size: 100 mL

lot number: \_\_\_\_\_

**Ingredients**

	final concentration	amount
RM059 sodium acetate, anhydrous	2.0 M	16.4 ± 0.4 g

**Procedure**

Slowly add the sodium acetate to approximately 50 mL deionized water.

Mix well.

Bring up to volume with deionized water.

Mix well.

Dispense into 100 mL bottles.

Autoclave at 250°F for 30 minutes.

Dispense into 15 mL tubes.

Store at room temperature.

**Data Log**

source	lot	amount
--------	-----	--------

RM059 sodium acetate, anhydrous	_____	_____
------------------------------------	-------	-------

**Quality Control**

QC023 QuantiBlot Quality Control of Solutions- Test 20 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCJ

Date: 10/1/96

**S044 Sodium Acetate, 0.2M**  
standard batch size: 250 mL

lot number: \_\_\_\_\_

**Ingredients**

	final concentration	amount
RM059 sodium acetate, anhydrous	0.2 M	4.1 ± 0.1 g

**Procedure**

Slowly add the sodium acetate to approximately 200 mL deionized water.

Mix well.

Bring up to volume with deionized water.

Mix well.

Dispense into 100 mL bottles.

Autoclave at 250°F for 30 minutes.

Store at room temperature.

**Data Log**

	source	lot	amount
RM059 sodium acetate, anhydrous	_____	_____	_____

**Quality Control**

QC023 QuantiBlot Quality Control of Solutions- Test 150 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RUDate: 10/27/96**S035 SP, 25X**

standard batch size: 4 L

lot number: \_\_\_\_\_

**Ingredients**

	final concentration	amount
RM003 EDTA	25 mM	37.2 ± 0.8 g
RM004 sodium hydroxide, 10N	-----	100 mL (guideline)
RM005 sodium chloride	3.75 M	877 ± 1 g
RM006 sodium phosphate, monobasic	0.25 M	138 ± 3 g

**Procedure**

Dissolve the EDTA in approximately 2 liters deionized water.

Adjust the pH to approximately 8.0 with 10N NaOH to help dissolve the EDTA.

Add the sodium phosphate. Mix until dissolved.

Add the sodium chloride.

Adjust the pH to  $7.4 \pm 0.2$  with 10N NaOH (about 80 mL).

Adjust the final volume to 4 liters with deionized water.

Measure and record the final pH.

Sterile filter each 2 liters with a sterile unexpired cellulose nitrate filter.

Store at room temperature.

**Data Log**

	source	lot	amount
RM003 EDTA	_____	_____	_____
RM004 sodium hydroxide, 10N	_____	_____	_____
RM005 sodium chloride	_____	_____	_____
RM006 sodium phosphate, monobasic	_____	_____	_____

**Quality Control**final pH: \_\_\_\_\_ specification  $7.4 \pm 0.2$ 

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *RG*

Date: *10/7/96*

**S108 SP, 2X**

lot number: \_\_\_\_\_

standard batch size: 1 L

**Ingredients**

final  
concentration

amount

S035 SP, 25X

2.0 X

80.0 ± 0.8 mL

**Procedure**

Add the SP to approximately 800 mL deionized water.

Bring up to the final volume with deionized water.

Dispense into a 500 mL bottles.

Autoclave at 250°F for 20 minutes.

Store at room temperature.

**Data Log**

source

lot

amount

S035 SP, 25X

\_\_\_\_\_

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCJ

Date: 10/2/96

**S098 Spotting Solution**

lot number: \_\_\_\_\_

standard batch size: 75 mL

Ingredients	final concentration	amount
S097 Pre-Wetting Solution	---	74.85 mL $\pm$ 1 mL
RM443 Bromothymol Blue, 0.04%	0.00008%	150 $\mu$ L $\pm$ 1 $\mu$ L

**Procedure**

Measure 74.85 mL Pre-Wetting Solution into a graduated cylinder and pour into a 100 mL bottle.

Add 150  $\mu$ L bromothymol blue.

Cap and mix well by inverting.

Store at room temperature.

**Data Log**

	source	lot	amount
S097 Pre-Wetting Solution	_____	_____	_____
RM443 Bromothymol Blue, 0.04%	_____	_____	_____

made by: \_\_\_\_\_

date: \_\_\_\_\_

Initials: 24

Date: 10/7/86

**S002 SSPE, 20X**

standard batch size: 4 L

lot number: \_\_\_\_\_

**Ingredients**

	final concentration	amount
RM003 EDTA	20. mM	29.8 ± 0.7 g
RM004 sodium hydroxide, 10N	-----	40 ± 5 mL (guideline)
RM005 sodium chloride	3.6 M	840 ± 10 g
RM006 sodium phosphate, monobasic	200 mM	110 ± 3 g

**Procedure**

Dissolve the EDTA in approximately 3 liters deionized water.

Adjust the pH to approximately 6.0 with 10N sodium hydroxide to help dissolve the EDTA.

Add the sodium phosphate first and then the sodium chloride.

Adjust the pH to 7.4 with 10N sodium hydroxide (about 40 mL).

Adjust the final volume to 4 liters with deionized water.

Measure and record the final pH.

Dispense into 1 L bottles.

Store at room temperature.

**Data Log**

	source	lot	amount
RM003 EDTA	_____	_____	_____
RM004 sodium hydroxide, 10N	_____	_____	_____
RM005 sodium chloride	_____	_____	_____
RM006 sodium phosphate, monobasic	_____	_____	_____

**Quality Control**

final pH: \_\_\_\_\_ specification 7.4 ± 0.2

made by: \_\_\_\_\_ date: \_\_\_\_\_



Initials: RU

Date: 10/7/96

# **S047 Stain Extraction Buffer**

lot number: \_\_\_\_\_

standard batch size: 1 L

## **Ingredients**

	final concentration	amount
S009 EDTA, 0.5M	10. mM	20 ± 1 mL
S052 TRIS-HCl, 0.1M - pH 7.8	10. mM	100 ± 0.5 mL
RM005 sodium chloride	100 mM	5.8 ± 0.2 g
RM101 dithiothreitol	33.9 mM	5.227 ± 0.008 g
S046 SDS, 20%	2.0 %	100 ± 3 mL
RM004 sodium hydroxide, 10N	-----	-----

## **Procedure**

Add all the ingredients except for the SDS to approximately 400 mL deionized water. Mix well.

Adjust the pH to 8.0 with 10N NaOH.

Record the pH.

Add the SDS. Mix well.

Bring up to the final volume with deionized water.

Dispense into sterile 125 mL bottles.

Store at 2-8°C.

## **Data Log**

	source	lot	amount
S009 EDTA, 0.5M	_____	_____	_____
S052 TRIS-HCl, 0.1M - pH 8.0	_____	_____	_____
RM005 sodium chloride	_____	_____	_____
RM101 dithiothreitol	_____	_____	_____
S046 SDS, 20%	_____	_____	_____
RM004 sodium hydroxide, 10N	_____	_____	_____

## **Quality Control**

final pH: \_\_\_\_\_ specification 8.0 ± 0.2

QC023 QuantiBlot Quality Control of Solutions- Test 150 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCJ

Date: 10/7/01

**S059 Sterile Water**

lot number: \_\_\_\_\_

standard batch size: 500 mL

**Procedure**

Filter sterilize 500 mL of deionized water.

Aliquot 10 mL each into 15 mL centrifuge tubes.

Autoclave at 250°F for 30 minutes.

Store at room temperature.

**Quality Control**

QC023 QuantiBlot Quality Control of Solutions- Test 150 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: ACS

Date: 10/7/86

**S023 Stripping Solution**

standard batch size: 2 L

lot number: \_\_\_\_\_

**Ingredients**

	final concentration	amount
RM102 formamide	55 %	1100 ± 10 mL
S035 SP, 25X	2.0 X	160 ± 4 mL
S046 SLS, 20%	1.0 %	100 + 6 mL

**NOTE:** S001 SDS, 20% can be substituted for 20% SLS in this solution.

**Procedure**

Add the SP and formamide to 1280 mL deionized water.

Mix well.

Add the SDS and mix gently.

Dispense into a 4 L brown bottle.

Store at 2-8°C.

**Data Log**

	source	lot	amount
RM102 formamide	_____	_____	_____
S035 SP, 25X	_____	_____	_____
S046 SLS, 20%	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: PCJ

Date: 10/21/96

**S017 TRIS-HCl, 1M - PH 8.0**

lot number: \_\_\_\_\_

standard batch size: 500 mL

Ingredients	final concentration	amount
RM073 TRIS	1.00 M	60.5 ± 0.1 g
RM096 hydrochloric acid	-----	-----

**Procedure**

Add the TRIS to approximately 400 mL deionized water.

Mix well.

Adjust the pH to 8.0 with concentrated hydrochloric acid.

Bring up to final volume with deionized water.

Measure and record the final pH.

Prepare a 1:100 dilution (10 mM TRIS-HCl) by mixing 1 mL TRIS-HCl solution and 99 mL deionized water.

Measure and record the pH of the dilution.

Autoclave at 250°F for 20 minutes.

Store at room temperature.

Data Log	source	lot	amount
RM073 TRIS	_____	_____	_____
RM096 hydrochloric acid	_____	_____	_____

**Quality Control**

final pH: \_\_\_\_\_ spec: 8.0 ± 0.1      1:100 pH: \_\_\_\_\_ spec: 8.0 ± 0.1

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RG

Date: 10/7/96

**S039 TE, 1X**

standard batch size: 500 mL

lot number: \_\_\_\_\_

**Ingredients**

	final concentration	amount
S049 TE, 100X	1.0 X	5.0 ± 0.3 mL

**Procedure**

Add the TE to approximately 400 mL deionized water.

Bring up to the final volume with deionized water.

Dispense into 125 mL bottles.

Autoclave at 250°F for 20 minutes.

Store at room temperature.

**Data Log**

	source	lot	amount
S049 TE, 100X	_____	_____	_____

**Quality Control**

final pH: \_\_\_\_\_ specification: 8.0 ± 0.2

QC023 QuantiBlot Quality Control of Solutions- Test 150 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCJ

Date: 10/21/96

**S128 TE<sup>-4</sup>, 1X**

lot number: \_\_\_\_\_

standard batch size: 500 mL

**Ingredients**

final  
concentration

amount

S017 TRIS-HCl, pH 8.0, 1 M

10 mM

5.0 ± 0.3 mL

S009 EDTA, 0.5 M

0.1 mM

100 ± 2 µL

**Procedure**

Add the TRIS and EDTA to 495 mL deionized water. Mix well

Filter Sterilize

Dispense into 15 mL sterile centrifuge tubes.

Autoclave at 250°F for 20 minutes.

Store at room temperature.

**Data Log**

source

lot

amount

S017 TRIS-HCl, pH 8.0, 1 M

\_\_\_\_\_

S009 EDTA, 0.5 M

\_\_\_\_\_

**Quality Control**

final pH: \_\_\_\_\_ specification: 8.0 ± 0.2

QC023 QuantiBlot Quality Control of Solutions- Test 150 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: Ad

Date: 10/7/96

**S049 TE, 100X**

standard batch size: 250 mL

lot number: \_\_\_\_\_

**Ingredients**

	final concentration	amount
RM003 EDTA	0.10 M	9.3 ± 0.5 g
RM073 TRIS	1.00 M	30.3 ± 0.1 g
RM004 sodium hydroxide, 10N	-----	-----
RM096 hydrochloric acid	-----	-----

**Procedure**

Add the EDTA to approximately 200 mL deionized water.

Adjust the pH to approximately 8.0 with sodium hydroxide to get the EDTA into solution. Mix until totally dissolved.

Add the TRIS and mix well.

Use hydrochloric acid or sodium hydroxide to adjust the pH of the solution to 8.0.

Bring up to final volume with deionized water.

Measure and record the final pH.

Dispense into 125 mL bottles.

Autoclave at 250°F for 30 minutes.

Store at room temperature.

**Data Log**

	source	lot	amount
RM003 EDTA	_____	_____	_____
RM073 TRIS	_____	_____	_____
RM004 sodium hydroxide, 10N	_____	_____	_____
RM096 hydrochloric acid	_____	_____	_____

**Quality Control**

final pH: \_\_\_\_\_ specification: 8.0 ± 0.2

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RS

Date: 10/7/96

# **S050 Test Gel Loading Buffer**

standard batch size: 100 mL

lot number: \_\_\_\_\_

## **Ingredients**

	final concentration	amount
RM020 bromophenol blue	0.10%	0.10 ± 0.01 g
RM217 xylene cyanol	0.10%	0.10 ± 0.01 g
RM040 Ficoll 400	5.0%	5.0 ± 0.1 g
S009 EDTA, 0.5M	20. mM	2.00 ± 0.05 mL
RM083 TAE, 10X	2.0 X	20.0 ± 0.5 mL

## **Procedure**

Combine the TAE, EDTA, and Ficoll.

Mix well. The solution may need to be heated gently to dissolve the Ficoll.

Add the bromophenol blue and xylene cyanol.

Mix well.

When all the solids are dissolved, bring up to volume using deionized water.

Filter sterilize.

Dispense 1.5 mL aliquots into 1.5 mL eppendorf tubes.

Store at -20°C.

## **Data Log**

	source	lot	amount
RM020 bromophenol blue	_____	_____	_____
RM217 xylene cyanol	_____	_____	_____
RM040 Ficoll 400	_____	_____	_____
S009 EDTA, 0.5M	_____	_____	_____
RM083 TAE, 10X	_____	_____	_____

made by: \_\_\_\_\_ date: \_\_\_\_\_



Initials: *RCJ*Date: *10/7/85***S062 Test Gel Standard**

lot number: \_\_\_\_\_

page 1 of 2

INGREDIENTS	initial concentration (ng/μL)	initial volume (μL)	final concentration	final volume (μL)
RM242 K562 DNA, Hae III fragments			5 ng/μL	
S050 test gel loading buffer	2 X		1 X	----
S059 sterile water	-----		----	----

**Calculations**

Record the initial concentration in ng/μL and the initial volume in μL of the K562 DNA, Hae III fragments received from the manufacturer.

Calculate the final volume according to equation 1.

$$(\text{final volume}) = \frac{(\text{initial DNA concentration})(\text{initial DNA volume})}{(5 \text{ ng/}\mu\text{L})} \quad \text{equation 1}$$

Record the final volume above. The final volume is the total batch size.

Calculate the amount of buffer to be added according to equation 2.

$$(\text{buffer volume}) = 0.5(\text{final volume}) \quad \text{equation 2}$$

Calculate the amount of sterile water to be added according to equation 3.

$$(\text{water volume}) = [0.5 * (\text{final volume})] - (\text{initial DNA volume}) \quad \text{equation 3}$$

Record the buffer and water volumes above.

To check the calculations, add together the initial volumes of DNA, loading buffer, and sterile water.

The sum of the initial volumes must be equal to the calculated final volume.

Initials: RU

Date: 10/17/96

**S062 Test Gel Standard**

lot number: \_\_\_\_\_

page 2 of 2

**Procedure**

Combine the DNA, loading buffer, and sterile water.

Mix well.

Using sterile pipet tips, dispense 500  $\mu$ L aliquots into sterile 1.5 mL eppendorf tubes.

Store at -20°C.

**Data Log**

	source	lot	amount
--	--------	-----	--------

RM242 K562 DNA, Hae III fragments	_____	_____	_____
--------------------------------------	-------	-------	-------

S050 test gel loading buffer	_____	_____	_____
---------------------------------	-------	-------	-------

S059 sterile water	_____	_____	_____
--------------------	-------	-------	-------

**Quality Control**

QC026 Gel Electrophoresis

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: PCSDate: 10/2/96**S051 TNE, 10X**

lot number: \_\_\_\_\_

standard batch size: 100 mL

**Ingredients**

	final concentration	amount
RM073 TRIS	100 mM	1.2 ± 0.02 g
S009 EDTA, 0.5M	10 mM	2.0 ± 0.1 mL
RM005 sodium chloride	1.0 M	5.8 ± 0.2 g
RM096 hydrochloric acid	---	

**Procedure**

Add the TRIS, EDTA, and sodium chloride to approximately 75 mL deionized water.  
Mix well.

Adjust the pH to 7.4 with hydrochloric acid

Bring up to the final volume with deionized water.

Measure and record the final pH. Adjust with concentrated HCl if necessary.

Filter to remove any particulates.

Dispense into a sterile 125 mL bottle.

Store at room temperature.

**Data Log**

	source	lot	amount
--	--------	-----	--------

RM073 TRIS

_____	_____	_____
-------	-------	-------

S009 EDTA, 0.5M

_____	_____	_____
-------	-------	-------

RM005 sodium chloride

_____	_____	_____
-------	-------	-------

RM096 hydrochloric acid

_____	_____	_____
-------	-------	-------

**Quality Control**

final pH: \_\_\_\_\_ specification: 7.4 ± 0.1

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: 24

Date: 10/27/96

**S107 TNE, 1X**

lot number: \_\_\_\_\_

standard batch size: 100 mL

**Ingredients**

final  
concentration

amount

S051 TNE, 10X

1.0 X

10.0 ± 0.3 mL

**Procedure**

Add the TNE to approximately 80 mL deionized water.

Bring up to the final volume with deionized water.

Dispense into a 125 mL bottles.

Autoclave at 250°F for 20 minutes.

Store at room temperature.

**Data Log**

source

lot

amount

S051 TNE, 10X

\_\_\_\_\_

**Quality Control**

QC023 QuantiBlot Quality Control of Solutions- Test 150 µL of solution

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCJ

Date: 10/7/96

**S052 TRIS, 0.1M - pH 7.8**

standard batch size: 1 L

lot number: \_\_\_\_\_

**Ingredients**

	final concentration	amount
RM073 TRIS	0.1 M	12.1 ± 0.2 g
RM096 hydrochloric acid	---	

**Procedure**

Add the TRIS to approximately 750 mL deionized water.

Mix well.

Adjust the pH to 7.8 with hydrochloric acid.

Bring up to the final volume with deionized water.

Mix well.

Dispense into a 1 L bottles.

Store at room temperature.

**Data Log**

	source	lot	amount
RM073 TRIS	_____	_____	_____
RM096 hydrochloric acid	_____	_____	_____

**Quality Control**

final pH: \_\_\_\_\_ spec: 7.8 ± 0.1

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RG

Date: 10/7/96

**S114 Urea (25 g Aliquot)**

lot number: \_\_\_\_\_

standard batch size: ~ 25 tubes x 25 g

**Ingredients**

aliquot

total amount

RM0079 Urea

25  $\pm$  0.2 g

625  $\pm$  6 g

electrophoresis grade

**Procedure**

**NOTE: WHEN WORKING WITH POWDERED UREA WEAR GLOVES, EYE PROTECTION, LAB COAT, AND RESPIRATOR FOR SAFETY.**

Fill out chemical logbook.

Using small weigh boat, weigh 25  $\pm$  0.2 g aliquots of urea.

Transfer the aliquots to 50 mL conical tubes

Cap all tubes tightly and label rack containing tubes with contents, lot number, date, initials, and safety data.

Store at room temperature.

**Data Log**

source

lot

amount

RM0079 urea

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: PCJ

Date: 2/9/98

**S114a Urea (18 g Aliquot-377 Sequencer)**

lot number: \_\_\_\_\_

standard batch size: ~ 25 tubes x 18 g

**Ingredients**

aliquot

total amount

RM0079 Urea

$18 \pm 0.1$  g

$450 \pm 4$  g

electrophoresis grade

**Procedure**

**NOTE: WHEN WORKING WITH POWDERED UREA WEAR GLOVES, EYE PROTECTION, LAB COAT, AND RESPIRATOR FOR SAFETY.**

Fill out chemical logbook.

Using small weigh boat, weigh  $18 \pm 0.1$  g aliquots of urea.

Transfer the aliquots to 50 mL conical tubes

Cap all tubes tightly and label rack containing tubes with contents, lot number, date, initials, and safety data.

Store at room temperature.

**Data Log**

source

lot

amount

RM0079 urea

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RCJ

Date: 2/9/98

**S114b Urea (10.8 g Aliquot-377 Sequencer)**

lot number: \_\_\_\_\_

standard batch size: ~ 25 tubes x 10.8 g

**Ingredients**

aliquot

total amount

RM0079 Urea

$10.8 \pm 0.1$  g

$450 \pm 4$  g

electrophoresis grade

**Procedure**

**NOTE: WHEN WORKING WITH POWDERED UREA WEAR GLOVES, EYE PROTECTION, LAB COAT, AND RESPIRATOR FOR SAFETY.**

Fill out chemical logbook.

Using small weigh boat, weigh  $10.8 \pm 0.1$  g aliquots of urea.

Transfer the aliquots to labeled 50 mL conical tubes.

Cap all tubes tightly and label rack containing tubes with contents, lot number, date, initials, and safety data.

Store at room temperature.

**Data Log**

source

lot

amount

RM0079 urea

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

made by: \_\_\_\_\_ date: \_\_\_\_\_



Initials: RCJ

Date: 10/7/96

**S024 Wash Solution #1**

lot number: \_\_\_\_\_

standard batch size: 4 L

**Ingredients**

final  
concentration

amount

S035 SP, 25X

2.0 X

320 ± 8 mL

S046 SLS, 20%

0.50 %

100 ± 1 mL

**NOTE:** S001 SDS, 20% can be substituted for 20% SLS in this solution.

**Procedure**

Add the SP to approximately 3 L deionized water.

Add the SLS to the solution.

Mix gently.

Bring up to the final volume with deionized water.

Dispense into a 4 L bottle.

Store at room temperature.

**Data Log**

source

lot

amount

S035 SP, 25X

\_\_\_\_\_

S046 SLS, 20%

\_\_\_\_\_

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: 20

Date: 10/7/96

**S025 Wash Solution #2**

lot number: \_\_\_\_\_

standard batch size: 4 L

**Ingredients**

final  
concentration

amount

S035 SP, 25X

0.10 X

16.0 ± 0.8 mL

S046 SLS, 20%

0.50 %

100 ± 1 mL

**NOTE:** S001 SDS, 20% can be substituted for 20% SLS in this solution.

**Procedure**

Add the SP to approximately 3 L deionized water.

Add the SLS to the solution.

Mix gently.

Bring up to the final volume with deionized water.

Dispense into a 4 L bottle.

Store at room temperature.

**Data Log**

source

lot

amount

S035 SP, 25X

\_\_\_\_\_

S046 SLS, 20%

\_\_\_\_\_

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: *RCW*Date: *10/21/96***S020 Yield Calibrators**

lot number: \_\_\_\_\_

standard batch size: 5 X 400 $\mu$ L each

page 1 of 2

**Ingredients**

	final concentration	amount
S128 TE <sup>-4</sup> , 1X	1 X	-----
RM148 lambda DNA	-----	140 $\pm$ 10 $\mu$ g (guideline)
S021 yield gel loading buffer	1.25 X	3.0 $\pm$ 0.5 mL
S059 Sterile Water	----	

**Calculations****Stock Solution**

Final DNA Concentration	Final Volume	Initial DNA Concentration	Volume Lambda DNA	Volume 1X TE <sup>-4</sup>
50 ng/ $\mu$ L	2800 $\mu$ L			

**Calibrators**

Calibrator	Final DNA Concentration	Stock DNA Concentration	Volume Stock DNA	Volume Water	Volume Buffer
A	300 ng/10 $\mu$ L	50 ng/ $\mu$ L	1200 $\mu$ L	300 $\mu$ L	500 $\mu$ L
B	200 ng/10 $\mu$ L	50 ng/ $\mu$ L	800 $\mu$ L	700 $\mu$ L	500 $\mu$ L
C	100 ng/10 $\mu$ L	50 ng/ $\mu$ L	400 $\mu$ L	1100 $\mu$ L	500 $\mu$ L
D	50 ng/10 $\mu$ L	50 ng/ $\mu$ L	200 $\mu$ L	1300 $\mu$ L	500 $\mu$ L
E	25 ng/10 $\mu$ L	50 ng/ $\mu$ L	100 $\mu$ L	1400 $\mu$ L	500 $\mu$ L
F	10 ng/10 $\mu$ L	50 ng/ $\mu$ L	40 $\mu$ L	1460 $\mu$ L	500 $\mu$ L

**Procedure**

Each lot of yield calibrators is prepared as a batch of five sets. Each batch requires 2800  $\mu$ L of 50 ng/ $\mu$ L stock lambda DNA solution.

Record the concentration in ng/ $\mu$ L of the lambda DNA received from the manufacturer under initial DNA concentration.

Initials: RO

Date: 10/21/96

## S020 Yield Calibrators

lot number: \_\_\_\_\_  
page 2 of 2

### Procedure

Calculate the volume of lambda DNA required for the stock solution according to equation 1.

$$(\text{volume lambda DNA}) = \frac{(\text{final DNA concentration})(\text{final volume})}{(\text{initial DNA concentration})} \quad \text{equation 1}$$

Calculate the volume of 1X TE<sup>-4</sup> to add to the stock solution according to equation 2.

$$(\text{volume 1X TE}^{-4}) = (\text{final volume}) - (\text{volume lambda DNA}) \quad \text{equation 2}$$

Prepare the stock solution by diluting the lambda DNA in a sterile centrifuge tube with 1X TE<sup>-4</sup> and mix well.

Label six sterile eppendorf tubes, one for each of the six yield calibrator levels.

Pipet the appropriate amounts of DNA stock solution and sterile water into the labeled tubes.

The combined volume of DNA and water is 1500 µL for each level.

Mix well.

Divide each level into five 300 µL aliquots, and dispense into labeled, sterile eppendorf tubes.

Add 100 µL of yield gel loading buffer to each tube. The final volume of each aliquot is 400 µL.

Store at -20°C.

### Data Log

S128 TE<sup>-4</sup>, 1X

RM148 lambda DNA

S021 yield gel loading buffer

S059 sterile water

source	lot	amount
_____	_____	_____
_____	_____	_____
_____	_____	_____
_____	_____	_____

### Quality Control

QC026 Gel Electrophoresis

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RC

Date: 10/7/86

# **S020 Yield Calibrators**

standard batch size: 5 X 400 $\mu$ L each

lot number: \_\_\_\_\_

page 1 of 2

## **Ingredients**

S039 TE, 1X  
RM148 lambda DNA  
S021 yield gel loading buffer  
S059 Sterile Water

final  
concentration  
1 X

amount

-----  
1.25 X

-----  
140  $\pm$  10  $\mu$ g (guideline)  
3.0  $\pm$  0.5 mL

## **Calculations**

Stock Solution				
Final DNA Concentration	Final Volume	Initial DNA Concentration	Volume Lambda DNA	Volume 1X TE
50 ng/ $\mu$ L	2800 $\mu$ L			

Calibrators					
Calibrator	Final DNA Concentration	Stock DNA Concentration	Volume Stock DNA	Volume Water	Volume Buffer
A	300 ng/10 $\mu$ L	50 ng/ $\mu$ L	1200 $\mu$ L	300 $\mu$ L	500 $\mu$ L
B	200 ng/10 $\mu$ L	50 ng/ $\mu$ L	800 $\mu$ L	700 $\mu$ L	500 $\mu$ L
C	100 ng/10 $\mu$ L	50 ng/ $\mu$ L	400 $\mu$ L	1100 $\mu$ L	500 $\mu$ L
D	50 ng/10 $\mu$ L	50 ng/ $\mu$ L	200 $\mu$ L	1300 $\mu$ L	500 $\mu$ L
E	25 ng/10 $\mu$ L	50 ng/ $\mu$ L	100 $\mu$ L	1400 $\mu$ L	500 $\mu$ L
F	10 ng/10 $\mu$ L	50 ng/ $\mu$ L	40 $\mu$ L	1460 $\mu$ L	500 $\mu$ L

## **Procedure**

Each lot of yield calibrators is prepared as a batch of five sets. Each batch requires 2800  $\mu$ L of 50 ng/ $\mu$ L stock lambda DNA solution.

Record the concentration in ng/ $\mu$ L of the lambda DNA received from the manufacturer under initial DNA concentration.

Initials: *ACJ*

Date: *10/17/86*

**S021 Yield Gel Loading Buffer**

lot number: \_\_\_\_\_

standard batch size: 100 mL

Page 1 of 2

**Ingredients**

final  
concentration

amount

RM020 bromophenol blue	0.25%	0.25 ± 0.01 g
RM217 xylene cyanol	0.25%	0.25 ± 0.01 g
RM040 Ficoll 400	12.5%	12.5 ± 0.1 g
S009 EDTA, 0.5M	50. mM	10.0 ± 0.1 mL
RM083 TAE, 10X	5.0 X	50.0 ± 0.5 mL
S001 SDS, 20%	0.20 %	1.00 ± 0.02 mL

**Procedure**

Combine the TAE, EDTA, SDS, and Ficoll.

Mix well. The solution may need to be heated gently to dissolve the Ficoll.

Add the bromophenol blue and xylene cyanol.

Mix well.

When all the solids are dissolved, bring up to volume using deionized water.

Filter sterilize.

Dispense 1.5 mL aliquots into sterile 1.5 mL eppendorf tubes.

Store at -20°C.

Initials: RY

Date: 10/7/95

**S021 Yield Gel Loading Buffer**

lot number: \_\_\_\_\_

standard batch size: 100 mL

Page 2 of 2

**Data Log**

source	lot	amount
--------	-----	--------

RM020 bromophenol blue	_____	_____
------------------------	-------	-------

RM217 xylene cyanol	_____	_____
---------------------	-------	-------

RM040 Ficoll 400	_____	_____
------------------	-------	-------

S009 EDTA, 0.5M	_____	_____
-----------------	-------	-------

RM083 TAE, 10X	_____	_____
----------------	-------	-------

S001 20% SDS	_____	_____
--------------	-------	-------

made by: \_\_\_\_\_ date: \_\_\_\_\_

Initials: RS

Date: 10/7/86

## S020 Yield Calibrators

lot number: \_\_\_\_\_  
page 2 of 2

### Procedure

Calculate the volume of lambda DNA required for the stock solution according to equation 1.

$$(\text{volume lambda DNA}) = \frac{(\text{final DNA concentration})(\text{final volume})}{(\text{initial DNA concentration})} \quad \text{equation 1}$$

Calculate the volume of 1X TE to add to the stock solution according to equation 2.

$$(\text{volume 1X TE}) = (\text{final volume}) - (\text{volume lambda DNA}) \quad \text{equation 2}$$

Prepare the stock solution by diluting the lambda DNA in a sterile centrifuge tube with 1X TE and mix well.

Label six sterile eppendorf tubes, one for each of the six yield calibrator levels.

Pipet the appropriate amounts of DNA stock solution and sterile water into the labeled tubes. The combined volume of DNA and water is 1500  $\mu$ L for each level.

Mix well.

Divide each level into five 300  $\mu$ L aliquots, and dispense into labeled, sterile eppendorf tubes.

Add 100  $\mu$ L of yield gel loading buffer to each tube. The final volume of each aliquot is 400  $\mu$ L.

Store at -20°C.

### Data Log

	source	lot	amount
S039 TE, 1X	_____	_____	_____
RM148 lambda DNA	_____	_____	_____
S021 yield gel loading buffer	_____	_____	_____
S059 sterile water	_____	_____	_____

### Quality Control

QC026 Gel Electrophoresis

made by: \_\_\_\_\_ date: \_\_\_\_\_