

FORENSIC BIOLOGY PROTOCOLS FOR FORENSIC STR ANALYSIS

Recovery of Spotted DNA Extracts from Whatman FTA Elute MicroCards

Status: Published

Document ID: 29291

DATE EFFECTIVE
09/20/2023

APPROVED BY
Nuclear DNA Technical Leader

PAGE
1 OF 3

Recovery of Spotted DNA Extracts from Whatman FTA Elute MicroCards

1 Procedure:

- 1.1 Obtain the spotted extract card(s) via the “Fbio_Spotted Extract Request” workflow in Qualtrax.
 - 1.1.1 The spotted card for the original extraction negative should also be pulled and reconstituted with the spotted sample.
- 1.2 Ad-Hoc samples into LIMS (See Creating Ad Hoc Samples in the LIMS Process Manual):
 - 1.2.1 If the samples were extracted and case handled pre-LIMS, create and shell in LIMS for the old case and ad-hoc the samples into LIMS as OC-Extracted DNA. Link the ad-hoc’ed samples to the old case shell.
- 1.3 Once all necessary samples are in LIMS, add them onto an “Extract-Recon” test batch for further processing.
- 1.4 The test batch should contain 1 new extraction negative for the recovery process (reconstitution extraction negative) and 1 for the spotted samples being reconstituted (if requested), unless that original extraction negative was already previously reconstituted. Print the input labels and label the appropriate number of 1.5 mL tubes.
- 1.5 Have a witness verify the samples, including the labels on the spotted extract cards against the test batch input samples.
- 1.6 Clean a 3mm hole puncher by punching a blank FTA Card 3 times, followed by wiping the puncher with a lint free wipe that has been prewet with 70% Ethanol.
- 1.7 Use the 3mm hole punch to punch out sample into a labeled 1.5 mL tube based on the chart below. The entire stain may be punched as needed, regardless of the concentration. The entire stain should be punched for spotted negative controls. All cards should be resealed in a kapak pouch and returned to storage, including those where the entire stain was punched.

<i>Original Quantitation Value</i>	<i>Recommended # of 3 mm punches for Recovery</i>
1 ng/uL and above	At least 8 punches
100 pg/uL to 999 pg/uL	At least 10 punches
Less than 100 pg/uL	Punch entire stain

- 1.7.1 A maximum of ten punches per 1.5mL tube is suggested for enhanced recovery. If more than ten punches are needed, separate the punches into multiple 1.5mL tubes.

Controlled versions of Department of Forensic Biology Manuals only exist in the Forensic Biology Qualtrax software. All printed versions are non-controlled copies.

© NYC OFFICE OF CHIEF MEDICAL EXAMINER

FORENSIC BIOLOGY PROTOCOLS FOR FORENSIC STR ANALYSIS

Recovery of Spotted DNA Extracts from Whatman FTA Elute MicroCards		
Status: Published		Document ID: 29291
DATE EFFECTIVE 09/20/2023	APPROVED BY Nuclear DNA Technical Leader	PAGE 2 OF 3

- 1.7.2 If multiple tubes are being used for the reconstitution of one sample, create an output sample in the LIMS batch for each reconstitution tube used. In the output sample data entry, indicate the number of punches in each tube.
- 1.7.3 The reconstitution extract from these samples can be combined during the microcon.
- 1.8 Heat the tubes with the punches at 80°C for 20 minutes to increase binding of the DNA to the card to prevent loss during washing.
- 1.9 Wash Step: Add 500 µL of sterile/ultrapure water to each tube, vortex for 5 seconds.
- 1.9.1 Note: Make sure the punches are moving around during the vortex step to ensure all punches are being washed. You may need to add more water if there are a large number of punches. Otherwise, there is the potential for PCR inhibitors to be left behind.
- 1.10 Briefly centrifuge the tubes, remove, and discard the supernatant.
- 1.11 Add sterile/ultrapure water to each tube, 50 µL per hole punch. Briefly centrifuge the tubes to ensure the punches are all immersed in water.
- 1.12 Place tubes in thermomixer at 95°C for 30 minutes, shaking at 500 rpm. At the end of the incubation step, remove the tubes from the heat block and vortex for 60 seconds.
- 1.13 Briefly centrifuge the tubes.
- 1.14 Collect the eluted DNA solution with a pipette and transfer to a 1.5 mL tube labeled with the LIMS output label. Discard the tube containing the discs/hole punches.
- 1.15 All reconstituted samples should then proceed to microcon via the Microcon testbatch. Samples should be microconned down to 25 µL. See [Microcon DNA Fast Flow DNA Concentration and Purification](#).
- 1.15.1 Samples can be microconned to concentrate, however, for samples where color from the card is left behind (usually samples with a greater number of punches), microcon to clean and concentrate.
- 1.15.2 Filter a maximum of 400µL at one time. For samples with more than 400µL, transfer the microcon membrane to a new collection tube as needed.
- 1.15.2.1 If samples were separated into multiple tubes during the reconstitution extraction, combine the reconstituted extracts using the same microcon membrane filter and continue to add the remaining extract to the microcon filter (up to 400µL each time). Additionally, under the input samples tab on the microcon batch, pool all the tubes from each sample

FORENSIC BIOLOGY PROTOCOLS FOR FORENSIC STR ANALYSIS

Recovery of Spotted DNA Extracts from Whatman FTA Elute MicroCards		
Status: Published		Document ID: 29291
DATE EFFECTIVE 09/20/2023	APPROVED BY Nuclear DNA Technical Leader	PAGE 3 OF 3

together, so there is one microcon output sample for each original sample card.

- 1.15.3 Store the final microconed extracts at 4°C and update the storage location in LIMS. The empty reconstitution tubes from the reconstitution test batch may be discarded.
- 1.15.4 After re-extraction, if sample testing will continue in Forensic Biology, a quantitation of the sample(s) is recommended.
- 1.15.5 If samples will be transferred outside of Forensic Biology, refer to [QC710 - Locating and Processing of Retained and/or Spotted Extracts](#)

ARCHIVED